**640200 : Potts Smith Dissecting Scissors**

Curved, 7 1/2" (19.0cm)

**354900 : DEBAKEY NEEDLE HOLDER**

13.7 x 2 mm jaw, 7" (18 cm), TC

**DP-60K : Teleflex® Premium Diamond Edge Aortic Punch - Standard Length****FREE-FLOAT CUTTING ACTION TECHNOLOGY** By allowing the diamond edge blade to "float" around the punch core during operation, alignment is maintained and mechanical jamming is prevented. One punch can be used to create...**P17258 : ANKENY STERNAL RETRACTOR**

6 sets of blades, 7" opening

**341290 : Davis-Finocchio Retractor, Infant**

Spread 2 3/4" (7cm), Side Blade 5/8" (1.6cm) wide x 1/2" (1.3cm) deep, Arm length 2 1/16" (5.2cm)

**354960 : Pilling® Thin Vascular and Coronary Technique® Needle Holder**

354960 Delicately proportioned, slender needle holders. 16.5 x 1.8 mm jaw, 6" (15 cm), TC




DP-56K : Teleflex® Premium Diamond Edge Aortic Punch - Standard Length

FREE-FLOAT CUTTING ACTION TECHNOLOGY By allowing the diamond edge blade to "float" around the punch core during operation, alignment is maintained and mechanical jamming is prevented. One punch can be used to create...



341250 : Burford-Finochietto Child Retractor

Spread 5 7/8" (15cm), Side Blade 15/16" (2.4cm), 1 5/16" (3.3cm), 1 3/4" (4.4cm), 1 3/4" (4.4cm) wide x 1" (2.5cm), 1 1/4" (3.2cm), 1 1/4" (3.2cm), 1 15/16" (4.9cm) deep, Arm length 4 5/16" (11cm)



354904 : DEBAKEY NEEDLE HOLDER

21.6 x 2.3mm jaw, 10 1/2" (26.5cm), TC



MDP-56K : Teleflex® Premium Diamond Edge Aortic Punch - Standard Length

FREE-FLOAT CUTTING ACTION TECHNOLOGY By allowing the diamond edge blade to "float" around the punch core during operation, alignment is maintained and mechanical jamming is prevented. One punch can be used to create...



341302 : Finochietto Retractor, Infant

Spread 2 3/4" (7cm), Side Blade 9/16" (1.5cm) wide x 1/2" (1.3cm) deep, Arm length 2 1/8" (5.5cm), Aluminum



342004 : WOLVEK STERNAL CRIMPER




841031 : KAISER-PILLING RING FORCEPS (2S)

Curved 3-3/8" jaw, 1 x 11 ratchets, 12-1/8" (31 cm)




341181 : Gerbode Sternal Retractor

Spread 8 1/2" (21.5cm), Side Blade 3 15/16" (10cm) wide x 1 3/16" (3cm) deep



341635 : Brewster Phrenic Retractor Nested Pair

341635 Nested set of 2 Retractors: Outside Retractor (double-ended) - overall length 149 mm (5-7/8") -Large End - 13 mm wide x 35 mm deep -Small End - 10 mm wide x 24 mm deep Inside Retractor (double-ended) - overall...




354966 : Micro Surgery Technique® Needle Holder

Flat handle and lock. Straight, 13 x 1.4 mm smooth jaw, 7" (18 cm), TC




342981 : DEBAKEY THORACIC TISSUE FORCEPS

Straight, 3mm tip, 7 3/4" (20 cm)



354916 : W.M.LEMMON NEEDLE HOLDER

1.5mm jaw, 6 3/4" (17cm), TC



357594 : Yasargil Vascular Scissors

Angled 45, bayonet, 7-1/2" (19 cm), 15 mm blades



841400 : K-P HOOK

Blunt, 13" (33 cm)



341773 : Harrington Splanchnic Retractor

Hollow handle, 2 1/2" (6.5cm) width x 5" (12.7cm) depth blade, 12" (30.5cm)



354974 : Castroviejo Needle Holder

Round handle and lock. Diama-grip jaws. Straight, jaw approximately 12 x 0.8 mm, 7-1/4" (18.5 cm)



341162 : Morse Sternal Retractor, Adult

Spread 8" (20cm), Side Blade 1 1/8" (2.8cm) wide x 15/16" (2.4cm) deep, Arm length 6" (15.2cm)



341225 : Burford Retractor, Adult

Stainless steel. Shanks are curved to fit body contour and blades are detachable. Blades have protruding points to help prevent spreader from moving during operation. Two pairs of blades are supplied. Minimum Spread:...



341330 : FINOCHIETTO RETRACTOR, MEDIUM

Spread 7 1/2" (19cm), Side Blade 2 5/8" (6.7cm) wide x 2 1/8" (5.4cm) deep, Curved arms