

**353641 : CRAWFORD ABDOMINAL CLAMP**

Straight, 6 3/8" (16.2cm) tip length, Angled shaft, 13 1/2" (34cm)

**353642 : CRAWFORD ABDOMINAL CLAMP**

Curved, 2 x 3 AtrauGrip Jaw, 3 1/4" (8.4cm) tip length, Angled shaft, 12 7/8" (32.5cm)

**353643 : CRAWFORD ABDOMINAL CLAMP**


Curved, 2 x 3 AtrauGrip Jaw, 2 5/8" (6.7cm) tip length, Curved shaft, 13 5/8" (34.5cm)

**353645 : Favalaro-Semb Ligature Carrier**


Jaws are constructed so that only the distal DeBakey- type Surfaces approximate. The smooth portion of the jaws remain open to prevent inadvertent clamping of the aortic wall. Half curved, DeBakey Teeth, 11-1/4" (28.5 cm)...

**353646 : Favalaro-Semb Ligature Carrier**

Jaws are constructed so that only the distal DeBakey- type Surfaces approximate. The smooth portion of the jaws remain open to prevent inadvertent clamping of the aortic wall. Full curved, DeBakey Teeth, 10-7/8" (27.5 cm)...


**353647 : DARLING TAPE PASSING CLAMP**

Curved, DeBakey Jaw, 2" (5cm) tip length, Curved shaft, 11 3/4" (29.5cm)



353650 : DEBAKEY AORTA CLAMP

2 x 3 AtrauGrip Jaw, 2 1/2" (6.5cm) tip length, Curved shaft, 10 1/4" (26cm)




353651 : GRANT AORTA CLAMP

2 x 3 AtrauGrip Jaw, 2 5/8" (6.7cm) tip length, Angled shaft, 7 7/8" (20cm)



353653 : GRANT AORTA CLAMP

Curved, 2 x 3 AtrauGrip Jaw, 2 1/8" (5.4cm) tip length, Angled shaft, 9 1/4" (23.5cm)




353655 : WYLIE J AORTIC CLAMP

Curved, DeBakey Jaw, 1 3/8" (3.65cm) tip length, Curved shaft, 9 5/8" (24.5cm)




353656 : WYLIE J AORTIC CLAMP

Curved, DeBakey Jaw, 1 5/8" (4.1cm) tip length, Curved shaft, 9 1/4" (23.5cm)




353657 : WYLIE HYPOGASTRIC CLAMP

Curved, DeBakey Jaw, 2.8mm tip width, Angled shaft, 9 5/8" (24.5cm)




353658 : WYLIE ABDOMINAL AORTA CLAMP

Curved, DeBakey Jaw, 1 1/8" (3cm) tip length, Angled shaft, 8 1/2" (21.5cm)



353660 : DEBAKEY TANGENTIAL CLAMP

Curved, 1 3/8" (3.5cm) tip length, Straight shaft, 7 3/4" (19.5cm)




353662 : DEBAKEY-TANGENTIAL CLAMP

Curved, 1 7/8" (4.76cm) tip length, Curved shaft, 10" (25.5cm)



353664 : DEBAKEY-TANGENTIAL CLAMP

Curved, 2" (5.1cm) tip length, Curved shaft, 10 3/8" (26.5cm)




353666 : DEBAKEY-TANGENTIAL CLAMP

Curved, 2 3/8" (5.9cm) tip length, Curved shaft, 10" (25.5cm)



353668 : DEBAKEY-TANGENTIAL CLAMP

Curved, 3 1/8" (7.8cm) tip length, Curved shaft, 10 1/8" (25.5cm)




353670 : DeBakey® Tangential Occlusion Clamps

Set of 5 instruments varying in length and jaw size. Dull, non-glare finished stainless steel. Includes: • 353660, No. 1, 1-3/8" (3.5 cm) jaw, 1 x 9 ratchets, 7-3/4" (19.5 cm) • 353662, No. 2, 1-7/8" (5 cm) jaw, 1 x 11...



353672 : STONEY AORTA CLAMP

Curved, DeBakey Jaw, 1/2" (1.4cm) tip length, Angled shaft, 6 7/8" (17.5cm)




353674 : STONEY AORTA CLAMP

Curved, DeBakey Jaw, 5/8" (1.45cm) tip length, Angled shaft, 9" (23cm)



353680 : DEBAKEY MULT-PURPOSE CLAMP

Angled 60, DeBakey Jaw, 1 1/4" (3.3cm) tip length, 9 1/2" (24cm)




353682 : DeBakey Multi-Purpose Angled Clamp

No. 2, angled 60, 1 x 9 ratchets, 8-5/8" (22 cm)



353684 : DeBakey Multi-Purpose Angled Clamp

Angled 60, 2" (5.08cm) tip length, 9 7/8" (25cm)



353686 : DeBakey Multi-Purpose Angled Clamp

Angled 60, 3 1/2" (8.9cm) tip length, 11 5/8" (29.5cm)