


CVT

**355834 : MARKING RING**


12 mm, 10/pk

**355838 : MARKING RING**


14 mm, 10/pk

**356025 : COOLEY PERIPHERAL VAS CLAMP**

Angled, 1 7/8" (4.9cm) tip length, 2.7mm tip width, Curved shaft, 7 1/2" (19cm)

**356027 : COOLEY ILIAC CLAMP**

Angled, 2 1/4" (5.7cm) tip length, Curved shaft, 8 1/4" (21cm)

**356028 : COOLEY ILIAC CLAMP**

Straight, 2 3/4" (7cm) tip length, Angled shaft, 9 3/4" (25cm)

**356031 : COOLEY PERIPHERAL VAS CLAMP**

Curved, 1 1/8" (2.9cm) tip length, 2mm tip width, Straight shaft, 6 3/4" (17cm)




356040 : COOLEY PAT-DUCTUS CLAMP

Straight, 2 1/8" (5.3cm) tip length, 2.3mm tip width, Straight shaft, 8" (20.5cm)




356041 : COOLEY PAT-DUCTUS CLAMP

Straight, 1 7/8" (4.9cm) tip length, 2.3mm tip width, Angled shaft, 7 7/8" (20cm)




356045 : COOLEY COARCTATION CLAMP

Straight, 2 7/8" (7.4cm) tip length, 2.4mm tip width, Straight shaft, 8 3/4" (22.5cm)




356046 : COOLEY COARCTATION CLAMP

Straight, 2 3/8" (6cm) tip length, 2.4mm tip width, Angled shaft, 8 7/8" (22.5cm)




356051 : COOLEY PARTIAL-OCCL CLAMP

Angled, 1" (2.6cm) tip length, 2.3mm tip width, Curved shaft, 8 1/8" (20.5cm)




356052 : COOLEY PARTIAL-OCCL CLAMP

Angled, 1 1/2" (3.9cm) tip length, 2.4mm tip width, Curved shaft, 8 7/8" (22.5cm)




356053 : COOLEY PARTIAL-OCCL CLAMP

Angled, 1 7/8" (4.8cm) tip length, 2.4mm tip width, Curved shaft, 8 7/8" (22.5cm)




356060 : COOLEY AORTA CLAMP

Strongly Curved, 2" (5.1cm) tip length, Straight shaft, 11 1/4" (28.5cm)




356070 : COOLEY MULTIPURPOSE CLAMP

Angled, 1 1/4" (3.3cm) tip length, 2.3mm tip width, Curved shaft, 9 1/8" (23cm)



356071 : COOLEY MULTIPURPOSE CLAMP

Angled, 1 5/8" (4.1cm) tip length, 2.3mm tip width, Curved shaft, 8 1/2" (21.5cm)



356072 : COOLEY MULTIPURPOSE CLAMP

Angled, 2 1/8" (5.4cm) tip length, 2.3mm tip width, Curved shaft, 8 1/2" (21.5cm)



357100 : PEDIATRIC CLAMP

Straight, Serrated, 1 1/2" (3.8cm) tip length, 1.1mm tip width, Straight shaft, 4 5/8" (12cm)



357102 : PEDIATRIC CLAMP

Angled, Serrated, 1 1/2" (3.8cm) tip length, 1.1mm tip width, Straight shaft, 4 5/8" (11.5cm)



357106 : PEDIATRIC CLAMP

Angled, Serrated, 1 1/8" (2.8cm) tip length, 1.1mm tip width, Straight shaft, 4" (10.5cm)



357110 : PEDIATRIC CLAMP

Slight Curved, 1/2" (1.18cm) tip length, 1.1mm tip width, Straight shaft, 4 1/2" (11.5cm)



357180 : CALMAN PEDIATRIC CLAMP

Angled, Serrated, 1" (2.6cm) tip length, 1.1mm tip width, Curved shaft, 4 1/2" (11.5cm)




357185 : EXAMPLE PEDIATRIC CLAMP

Angled, Serrated, 1" (2.5cm) tip length, 1.1mm tip width, Curved shaft, 5 5/8" (14cm)



357191 : CASTANEDA PEDIATRIC CLAMP

Curved, Serrated, 1/2" (1.4cm) tip length, 1.1mm tip width, Curved shaft, 6" (15cm)



357192 : CASTANEDA PEDIATRIC CLAMP

Curved, Serrated, 3/4" (2.03cm) tip length, 1.1mm tip width, Curved shaft, 6" (15.5cm)

