


## CVT

**357194 : THOMPSON PEDIATRIC CLAMP**

Angled, Serrated, 1/2" (1.15cm) tip length, 1.1mm tip width, Curved shaft, 6 1/8" (15.5cm)

**357195 : THOMPSON CAROTID CLAMP**

Angled, DeBakey Jaw, 1.6mm tip width, Curved shaft, 5 3/4" (14.5cm)

**357197 : KARTCHNER CAROTID CLAMP**

Angled, Serrated, 3/8" (1.1cm) tip length, 1.1mm tip width, Curved shaft, 6 3/8" (16.5cm)

**357198 : KARTCHNER CAROTID CLAMP**

Angled, DeBakey Jaw, 1.7mm tip width, Curved shaft, 6 1/4" (16cm)

**357200 : PEDIATRIC CLAMP**

Straight, Serrated, 1 7/8" (4.8cm) tip length, 1.1mm tip width, Straight shaft, 6 1/8" (15.5cm)

**357202 : PEDIATRIC CLAMP**

Angled, Serrated, 1/2" (1.27cm) tip length, 1.1mm tip width, Straight shaft, 6" (15.5cm)



**357262 : PEDIATRIC CLAMP**

Angled, Serrated, 3/4" (1.9cm) tip length, 1.1mm tip width, Curved shaft, 5 5/8" (14.5cm)



**357270 : PEDIATRIC CLAMP**

Angled, Serrated, 3/4" (1.8cm) tip length, 1.1mm tip width, Curved shaft, 5 7/8" (15cm)



**357302 : COOLEY NEONATAL CLAMP**

Straight, 2mm tip width, Angled shaft, 5 1/8" (13cm)



**357304 : COOLEY NEONATAL CLAMP**

Angled, 2mm tip width, Curved shaft, 5" (12.5cm)




**357306 : COOLEY NEONATAL CLAMP**

Curved, 3/4" (19mm) tip length, 2mm tip width, Straight shaft, 4 5/8" (11.5cm)



**357308 : COOLEY NEONATAL CLAMP**

Angled, 2mm tip width, Curved shaft, 4 5/8" (12cm)



**357310 : COOLEY NEONATAL CLAMP**

Angled, 2mm tip width, Curved shaft, 4 3/4" (12cm)

**357580 : MICRO BULLDOG CLAMP**

Straight, Cross-serrated, 1/4" (7mm) tip length, 1mm tip width, 3/4" (2cm)



**357581 : MICRO BULLDOG CLAMP**

Straight, Cross-serrated, 3/8" (8mm) tip length, 1.3mm tip width, 7/8" (2cm)




**357582 : MICRO BULLDOG CLAMP**

Angled, Cross-serrated, 3/8" (9mm) tip length, 1.4mm tip width, 7/8" (2cm)



**357583 : MICRO BULLDOG CLAMP**

Curved, Cross-serrated, 3/8" (9.5mm) tip length, 1.2mm tip width, 3/4" (2cm)



**357590 : MINI BULLDOG CLAMP**

Straight, Cross-serrated, 1/8" (4mm) tip length, 1mm tip width, 1/2" (1cm)



**357591 : MINI BULLDOG CLAMP**


Straight, Cross-serrated, 1/4" (6.3mm) tip length, 1mm tip width, 5/8" (1.5cm)



**357592 : MINI BULLDOG CLAMP**

Curved, Cross-serrated, 1/4" (7mm) tip length, 5/8" (1.5cm)





**357650 : Pilling Micro Vascular Set**



**357659 : MICRO TYING FORCEPS**

Curved, Tying Platform, Smooth, 0.5mm tip, 6" (15 cm)



**357661 : MICRO TYING FORCEPS**

Curved, Tying Platform, Smooth, 0.5mm tip, 7 1/8" (18 cm)



**357681 : MICRO GRASPING FORCEPS**

Straight, Serrated, 7 7/8" (20 cm)



**358007 : MICRO DENIS RING FORCEPS**

Straight, Smooth, 1mm tip, 8 1/4" (210 cm), Diama-Grip jaws