

CVT

**359021 : DEBAKEY TITANIUM TISSUE FORCEPS**

2.5mm tip, 6" (15 cm)

**359027 : DEBAKEY TITANIUM TISSUE FORCEPS**

Delicate, 3mm tip, 6" (15 cm)

**359028 : DEBAKEY THORACIC TISSUE FORCEPS**

3mm tip, 7 5/8" (20 cm)

**359030 : DEBAKEY TITANIUM TISSUE FORCEPS**

Angled, 2.5mm tip, 5 7/8" (15 cm)

**359032 : DEBAKEY TITANIUM TISSUE FORCEPS**

Angled, 2.5mm tip, 9 1/4" (24 cm)

**359033 : MICRO TITANIUM DIAMAGRIP TYING FORCEPS**

4 3/4" (12 cm), Diamo-Grip jaws



359034 : MICRO TITANIUM DIAMAGRIP TYING FORCEPS

4 3/4" (12 cm), Diama-Grip jaws



359036 : MICRO TITANIUM DIAMAGRIP TYING FORCEPS

Round handle. Curved, 5-7/8" (15 cm), 6 mm Diama-Grip jaws



359042 : MICRO TITANIUM DIAMAGRIP TYING FORCEPS

0.4mm tip, 7 1/8" (18 cm), Diama-Grip jaws



359043 : ADSON TITANIUM TISSUE FORCEPS

Serrated, 4 3/4" (12 cm), TC



359050 : MCINDOE TITANIUM TISSUE FORCEPS

Serrated, 5 7/8" (15 cm)



359051 : GILLES TITANIUM TISSUE FORCEPS

Serrated, 1 x 2 Teeth, 5 7/8" (15 cm)



359053 : BARRAYA TITANIUM TISSUE FORCEPS

Straight, Serrated, 2 x 3 Teeth, 7 3/4" (20 cm)

359054 : GUILMET TITANIUM DRESSING FORCEPS

Angled, Serrated, 7 7/8" (20 cm)



359122 : DEBAKEY TITANIUM BULLDOG CLAMP

Straight, 3/4" (2cm) tip length, 3" (7.5cm)



359123 : DEBAKEY TITANIUM BULLDOG CLAMP

Straight, 1 1/8" (3cm) tip length, 3 3/8" (8.5cm)



359124 : DEBAKEY TITANIUM BULLDOG CLAMP

Straight, 1 3/4" (4.4cm) tip length, 3 7/8" (10cm)



359127 : DEBAKEY TITANIUM BULLDOG CLAMP

Curved, 3 1/4" (8.5cm)



359128 : DEBAKEY TITANIUM BULLDOG CLAMP

Curved, 3 7/8" (10cm)



359131 : DEBAKEY TITANIUM BULLDOG CLAMP

Angled 45, 4 3/4" (12cm)





359138 : GLOVER TITANIUM COARCTATION CLAMP

Straight, DeBakey Jaw, 8 5/8" (22cm)




359139 : GLOVER TITANIUM COARCTATION CLAMP

Straight, DeBakey Jaw, Angled shaft, 8 5/8" (22cm)




359152 : KAY TITANIUM VASCULAR CLAMP

Curved, DeBakey Jaw, Angled shaft, 7 7/8" (20cm)



359159 : DEBAKEY TITANIUM VASCULAR CLAMP

Angled, 2 x 3 AtrauGrip Jaw, 2 1/2" (6.35cm) tip length, Curved shaft, 10 1/4" (26cm)



359160 : HUFNAGEL TITANIUM VASCULAR CLAMP

Angled, DeBakey Jaw, 3" (7.5cm) tip length, Angled shaft, 8 5/8" (22cm)