
**359215 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS
15 GRAMS**


Angled, DeBakey type jaw, 5 1/2" (14 cm)

**359221 : COOLEY TITANIUM PEDIATRIC CLAMP**


Spoon, Serrated, 2.2mm tip width, 7 1/4" (18.5cm)

**359236 : DERRA TITANIUM ANASTAMOSIS CLAMP**

Curved, DeBakey Jaw, 1 1/8" (2.7cm) tip length, 2.6mm tip width, 6 3/4" (17cm)

**359237 : DERRA TITANIUM ANASTAMOSIS CLAMP**

Curved, DeBakey Jaw, 2.8mm tip width, 6 3/4" (17.5cm)

**359238 : RENAL ARTERY TITANIUM CLAMP**

Angled, Serrated, 2.5mm tip width, 7 1/2" (19cm)

**359283 : MICRO TITANIUM DIAMAGRIP TYING FORCEPS**

Straight, I.D. 1.0mm tip, 8 5/8" (22 cm), Diama-Grip jaws



359296 : HALSTEAD TITANIUM MOSQUITO FORCEPS

Curved, Serrated, 4 7/8" (13 cm)




359307 : CRILE TITANIUM HEMOSTATIC FORCEPS

Curved, Serrated, 5 5/8" (15 cm)



359313 : O'SHAUGHNESSY TITANIUM FORCEPS

Right angle, Serrated, 8 1/4" (21 cm)



359314 : HEISS TITANIUM FORCEPS

Angled, Serrated, 8 1/4" (21 cm)




359321 : O'SHAUGHNESSY TITANIUM FORCEPS

Right angle, Serrated, 2.5mm tip, 9 1/2" (24 cm)



359325 : CRAFOORD TITANIUM CLAMP FORCEPS

Curved, Serrated, 9 1/2" (24 cm)



359354 : DEBAKEY TITANIUM PEDIATRIC CLAMP

Straight, 1.6mm tip width, 4 3/8" (11cm)

NO IMAGE
AVAILABLE

359355 : DEBAKEY TITNIUM PEDIATRIC CLAMP

Angled 60, 7/8" (2.3cm) tip length, 1.6mm tip width, Curved shaft, 4" (10.5cm)



359356 : DEBAKEY TITNIUM ANASTOMOSIS CLAMP

Angled, 3/8" (0.95cm) tip length, 1.6mm tip width, Angled shaft, 3 7/8" (10cm)



359357 : DEBAKEY TITNIUM ANASTOMOSIS CLAMP

Spoon, 1/4" (0.64cm) tip length, 1.2mm tip width, 3 7/8" (9.5cm)



359358 : DEBAKEY TITNIUM ANASTOMOSIS CLAMP

Angled, 5/8" (1.59cm) tip length, 1.6mm tip width, 4 1/8" (10.5cm)



359360 : CUSHING TITANIUM TISSUE FORCEPS

Serrated, 7 7/8" (20 cm)



359361 : ADSON TITANIUM DRESSING FORCEPS

Serrated, 4 3/4" (12 cm)




359383 : Finochietto Titanium Retractor

4 1/4" (10.5cm)




359407 : DENIS MICRO TITANIUM FORCEPS

Straight, 0.4mm tip, 8 1/4" (21 cm), Diama-Grip jaws



506743 : Mckenzie "V" Ligating Clips

Smooth, vial of 100, McKenzie Ligating Clips are non-magnetic and made of sterling silver



506744 : Mckenzie "V" Ligating Clips

Serrated, vial of 100, McKenzie Ligating Clips are non-magnetic and made of sterling silver



506746 : Duane "U" Ligating Clips

Smooth, vial of 100, Duane Ligating Clips are non-magnetic and made of sterling silver



506749 : Clip Holding Magazine For Mckenzie and Duane Clips

Single row, weighted, manual loading