

## CVT

NO IMAGE  
AVAILABLE

**P17951 : MICRO TISSUE FORCEPS**

Lightweight, w/ platform, 1 x 2, 7-1/4" (18.5 cm)

NO IMAGE  
AVAILABLE

**P18156 : DEBAKEY TITANIUM CLAMP**

Straight 3.5 cm jaw, angled handle

NO IMAGE  
AVAILABLE

**P18182 : COOLEY NEONATAL VASCULAR CLAMP**

Carotid clamp, 5" (12.5 cm)

NO IMAGE  
AVAILABLE

**P18183 : CLEVELAND CLINIC CLAMP, LEFT, LARGE**

9-5/8" (24.5 cm)

NO IMAGE  
AVAILABLE

**P18184 : CLEVELAND CLINIC CLAMP, RIGHT, LARGE**

9-5/8" (24.5 cm)

NO IMAGE  
AVAILABLE

**P18185 : KITZMILLER CLAMP**

Left angled jaw, 5-1/2" (14 cm)

NO IMAGE  
AVAILABLE

**P18186 : KITZMILLER CLAMP**

Right angled jaw, 5-1/2" (14 cm)

NO IMAGE  
AVAILABLE

**P18246 : CASTENADA CLAMP**

Curved shank, jaw angled 30 deg., 5" (12.5 cm)

NO IMAGE  
AVAILABLE

**P18250 : CASTENADA CLAMP, MEDIUM**

Curved shank, 6" (15 cm)

NO IMAGE  
AVAILABLE

**P18261 : SELMAN COOLEY FORCEPS, LARGE**

Angled, 5-1/2" (14 cm)

NO IMAGE  
AVAILABLE

**P18308 : COOLEY ILIAC CLAMP**

6 cm jaw, 10" (25.5 cm)

NO IMAGE  
AVAILABLE

**P18481 : FOGARTY HYDRAGRIP CLAMP**

Straight, 6-1/2" (16.5 cm)

NO IMAGE  
AVAILABLE

**P18482 : FOGARTY HYDRAGRIP CLAMP**

Angled jaws, 5" (12.5 cm)



**P18483 : FOGARTY HYDRAGRIP CLAMP**

Straight, 4" jaw, 11-1/2" (29 cm)



**P18484 : FOGARTY HYDRAGRIP CLAMP**

Angled jaws, 8-3/4" (22 cm)




**P18486 : HEMOGRIP CLAMP**

Angled, 3-3/4" jaw, 8-1/2" (21.5 cm)




**P18554 : WELDON MINI-BULLDOG CLAMP**

Straight, 1-3/8" (3.5 cm)




**P18555 : WELDON MINI-BULLDOG CLAMP**

Curved, 1-1/4" (3 cm)



**P18587 : CLAMP BOX**

70 x 45mm, stainless steel



**P18663 : CHEST TUBE PASSER**

9-3/4" (25 cm)

AVAILABLE

NO IMAGE  
AVAILABLE

**P18664 : CLAMP CHEST TUBE PASSER**

W/ratchet, 9-3/4" (25 cm)

NO IMAGE  
AVAILABLE


**P18665 : CLAMP CHEST TUBE PASSER**

14" (35.5 cm)



**P18667 : CRAWFORD-COOLEY TUNNELER**

12mm, 1.5" (4cm)



**P18675 : HALLMAN TUNNELING FORCEPS**

14" (35.5 cm)

NO IMAGE  
AVAILABLE

**P18697 : COOLEY CAVAL OCCLUSION FORCEPS**

36 fr., 8-1/4" (21 cm)