

## CVT



NO IMAGE  
AVAILABLE

**P18741 : PULMONARY ARTERY CLAMP**

12 mm, malleable shaft, 9" (23 cm)



NO IMAGE  
AVAILABLE

**P18743 : PULMONARY ARTERY CLAMP**


20 mm, malleable shaft, 9" (23 cm)



NO IMAGE  
AVAILABLE

**P18745 : PULMONARY ARTERY CLAMP**


30 mm, malleable shaft, 9" (23 cm)



NO IMAGE  
AVAILABLE

**P18777 : MYOCARDIAL DILATOR**

2mm tip graduated to 17.5mm, 5 3/4" (14.5cm)



NO IMAGE  
AVAILABLE

**P18871 : RUMEL TOURNIQUET FLEX-NECK, ADULT**

12" (30.5cm)



NO IMAGE  
AVAILABLE

**P18873 : RUMEL TOURNIQUET FLEX-NECK, ADULT**

w/ratchet, 12" (30.5cm)

NO IMAGE  
AVAILABLE

**P34409 : KITZMILLER TITANIUM RASP, LEFT ANGLED**

4 1/4" (11cm)

NO IMAGE  
AVAILABLE

**P34410 : KITZMILLER TITANIUM RASP, RIGHT ANGLED**

4 1/4" (11cm)

NO IMAGE  
AVAILABLE

**P35323 : In-Situ Titanium Bypass Hook**

Angled, ball tip, 7" (18cm)