

**841100 : KAISER-PILLING BIOPSY FORCEPS**

Malleable shaft, 4 mm cup, 12-5/8" (32 cm)

**341525 : Malleable Atrial Retractor**


3/4" (2.0cm) width blade, 8 3/4" (22.0cm)

**351909 : GERALD ATRAUGRIP TISSUE FORCEPS**


DeBakey type jaw, 9 1/2" (24 cm)

**351668 : RETRACTOR SIDE BLADE FOR 351675 PILLING MULTIPURPOSE RETRACTOR**

4 Blunt prongs, Side Blade 7/8" (2.2cm) wide x 3 1/8" (8.0cm) deep

**354906 : DEBAKEY NEEDLE HOLDER**

19 x 2.6mm jaw, 11 3/4" (30cm), TC


**355275 : HIEBERT DILATOR**

0.5mm x 1mm, Malleable shafts, 6" (15.5cm)

NO IMAGE
AVAILABLE

P17773 : VARADY EXTRACTOR

Double-ended, 7" (18 cm). Crochet hook on one end.



352165 : Pilling Coronary Potts Scissors

Angled 45, 7" (18.0cm), Delicate, 16.0mm blades



355091 : Garrett Vascular Dilator

Long 1.0 mm, malleable shaft *Related: 355080 Short 1.0 mm*
355081 Short 1.5 mm ◆ *Long 1.5 mm 355082 Short 2.0...*

◆ Long 1.0 mm



841013 : KAISER-PILLING RING FORCEPS (3L)

Long jaw, curved on side, 12-1/8" (31 cm)

NO IMAGE
AVAILABLE


P17981 : JACOBSON NEEDLE HOLDER

Straight, delicate, flat handle, w/lock, 8 1/4" (21cm)



211245A : LANDRENEAU LAWRENCE FINE DISSECTING FORCEPS

45 degree curve, locking, 21 cm shaft length, 15-3/8" (39 cm) overall length. Not designed for use with a 10 mm cannula.



355559DG : CASTROVIEJO MICRO NEEDLE HOLDER

Round handle, w/Lock, 1mm Diama-Grip jaw, 10" (25.5cm)



341389 : TRANSPLANT DONOR STERNAL RETRACTOR, ADULT

Spread 7 1/2" (19cm), Sharp, Side Blade 7/8" (2.3cm) wide x 1" (2.5cm) deep



341511 : LEMOLE-PILLING T FRAME ONLY FOR TYPE #1 RETRACTOR

Lemole-Pilling T-Frame #1 only



351804 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 9 1/2" (24 cm)




355092 : Garrett Vascular Dilator

Long 1.5 mm, malleable shaft *Related: 355080 Short 1.0 mm* ♦ Long 1.0 mm
355081 Short 1.5 mm ♦ Long 1.5 mm 355082 Short 2.0...




341650 : Davidson Scapula Retractor, Small

2" (5cm) width blade, 7 3/4" (19.5cm)



355278 : HIEBERT DILATOR

3.5mm x 4mm, Malleable shafts, 6" (15.5cm)



342557 : KANTROWITZ FORCEPS

Right Angle, Serrated, 3/4" (1.75cm) tip length, 2mm tip width, 7 7/8" (20cm)



341170 : Gerbode-Burford Retractor, Adult

Adult Spread 7 1/2" (19 cm), Small blade 2 1/8" deep x 2 5/8" wide (55 x 67 mm), Large blade 2 7/8" deep x 2 5/8" wide (73 x 67 mm), Arm length 7 1/4" (18.5 cm), Stainless steel




351800 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 6" (16 cm)



355302 : Debakey® Vascular Dilator

6.0 mm, 7-1/2" (19 cm) long *Related: 355295 1.0 mm, 7-1/2" (19 cm) long 355290 Set of 12, 15" (38 cm), in case 355296 1.5 mm, 7-1/2" (19 cm) long...*



358001 : JACOBSON MICRO NEEDLE HOLDER

Straight, Round handle, With lock, 11 x 0.8mm Diama-Grip jaw, 9" (23cm)



341303 : Amato Neonatal Sternal Retractor, small

Specially designed infant sternal retractors to reduce post-operative sternal deformity by dispensing pressure over greater surface area. Selection of 3 sizes available for proper fit. 4 blade, small Swivel blades 3/8"...

