


CVT

**354972 : Castroviejo Needle Holder**

Round handle and lock. Diama-grip jaws. Straight, jaw 11.7 x 0.8 mm, 6" (15 cm)

**355276 : HIEBERT DILATOR**


1.5mm x 2mm, Malleable shafts, 6" (15.5cm)

**357595 : Yasargil Vascular Scissors**

Reverse angled 125, Bayonet, 7 1/4" (18.5cm), 12.0mm blades

**841030 : KAISER-PILLING RING FORCEPS (1S)**

Straight, 3-3/8" jaw, 1 x 11 ratchets, 12-5/8" (32 cm)

**342572 : SAROT THORACIC ARTERY FORCEPS**

Curved, Serrated, 2 1/4" (5.7cm) tip length, 2mm tip width, 9" (23cm)

**341522 : Cooley Atrial Retractor, Small**

5/8" (1.5cm) width blade, 9 1/2" (24.0cm)



351672 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR

Sharp, Center blade 5/8" (1.6cm) wide x 2" (5.2cm) deep



352150 : Potts Smith Scissors

Angled 25, 7 1/2" (19.0cm), 25.0mm blades



351785 : DEBAKEY VALVE HOOKS

Small, straight, 3/8" (8.6 mm) hook, 12-1/2" (31.6 cm)




355080 : Garrett Vascular Dilator

Short 1.0 mm, malleable shaft *Related: 355080 Short 1.0 mm 355081 Short 1.5 mm*

◆ Long 1.0 mm

◆ Long 1.5 mm 355082 Short 2.0...



352154 : Potts Smith Scissors

Angled 60, 7" (18.0cm), 23.0mm blades



355198 : DOERING CAROTID END POINT RONGEUR FORCEPS

2.5mm cup jaws, 7 1/4" (19 cm)



341790 : Harrington Extra Deep Retractor

Hollow handle, 2 1/4" (5.7cm) width x 8" (20.5cm) depth blade, 11 3/4" (30.0cm)

342596 : GEMINI FORCEPS

Curved, Serrated, 1/2" (1.23cm) tip length, 1.3mm tip width, 7 1/4" (18.5cm)



P18875 : RUMEL TOURNIQUET FLEX-NECK OBTURATOR ONLY, ADULT

12" (30.5cm)



355277 : HIEBERT DILATOR

2.5mm x 3mm, Malleable shafts, 6" (15.5cm)



341506 : LEMOLE-PILLING® VALVE RETRACTOR - Set #1 with Cooley Retractor (341195)

Provides self-retaining valve retraction. Fits directly onto arms of many types of sternal retractors. The T-Frame cross bar attaches to the sternal retractor by means of a "C" clamp. 8" (20.3 cm) T-Frame cross bar allows...




351791 : Cooley Mitral Valve Retractor

Curved left, 2" (5cm) width blade, 8 7/8" (22.5cm)



352152 : Potts Smith Scissors

Angled 45, 7 1/4" (18.5cm), 23.0mm blades




355081 : Garrett Vascular Dilator

Short 1.5 mm, malleable shaft *Related: 355080 Short 1.0 mm*
355081 Short 1.5 mm

◆ Long 1.0 mm

◆ Long 1.5 mm 355082 Short 2.0...





359335 : MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, With lock, 9 x 1.2mm jaw, 7 1/8" (18cm), TC



354320 : PILLING RATCHET LOCK

5/8" (1.5cm)



359352 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, With lock, 1.5mm jaw, 9" (23cm), TC



841052 : KAISER-PILLING ATRAUGRIP RING FORCEPS (2AG)

Curved, DeBakey type jaw, 12-1/4" (31 cm)



841401 : K-P HOOK

Insulated, 13" (33 cm)