

NO IMAGE
AVAILABLE

518320 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #4, SHORT LENGTH

5.5mm I.D., 8mm O.D., 52mm, stainless steel

NO IMAGE
AVAILABLE

518245 : JACKSON IMPROVED TRACHEOSTOMY TUBE #9, SHORT LENGTH

9.9mm I.D., 13mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

507375 : JACKSON CROSS ACTION FORCEPS

Angled, 9-1/4" (24 cm)

NO IMAGE
AVAILABLE

516400 : MARTIN LARYNGECTOMY TUBE #12

12.7mm I.D., 16.3mm O.D., 42mm, sterling silver

NO IMAGE
AVAILABLE

P26467 : FOOTPLATE PICK

Angled up 90 deg., 6 3/8" (16cm)

NO IMAGE
AVAILABLE

073502 : Hurd Pillar Retractor, Large

10mm / 14mm width x 5mm depth blade, 9" (23.0cm)

NO IMAGE
AVAILABLE

506470 : Jako-Kleinsasser Micro-Laryngeal Grasping Forceps

These Pilling Micro-laryngeal Grasping Forceps were designed specifically for everyday-use micro-surgery. They feature slender 2 mm diameter stems with miniature distal jaws in the shape of grasping forceps and are made of...

NO IMAGE
AVAILABLE

507828 : JACKSON ESOPHAGEAL DILATOR

28 Fr.

NO IMAGE
AVAILABLE

518050 : Jackson Improved Tracheostomy Tube #9

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518050 (size 9, 13.0 mm O.D., 9.9 mm...

NO IMAGE
AVAILABLE

P26468 : FOOTPLATE PICK

Angled down 90 deg., 6 3/8" (16cm)

NO IMAGE
AVAILABLE

507810 : JACKSON ESOPHAGEAL DILATOR

10 Fr.

NO IMAGE
AVAILABLE

518330 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #6, SHORT LENGTH

7.2mm I.D., 10mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

511134P : TUCKER VALVE INNER TUBE ONLY FOR #5 JACKSON IMPROVED TUBE

6.1mm I.D., 71mm, sterling silver



514376 : JACKSON CHILD DILATOR

5 1/2" (14cm)



411165 : FARRELL APPLICATOR TWISTED

6 1/4" (16.0cm)



519005 : JACKSON LARYNGECTOMY TUBE #7

9.1mm I.D., 12mm O.D., 55mm, stainless steel



506862 : FEDER-OSSOFF MICRO LARYNGEAL CUP FORCEPS

Straight. Shaft tapers from 2.0 mm to 1 mm. Round biopsy cups w/ 1 mm jaw, 9" (23 cm)



518340 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #8, SHORT LENGTH

9.2mm I.D., 12mm O.D., 56mm, stainless steel



518100 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #00

2.4mm I.D., 4.5mm O.D., 33mm, stainless steel



507501 : JACKSON LARYNGEAL DILATOR

12 Fr., 8 1/2" (21.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

518044 : Jackson Improved Tracheostomy Tube - Fenestrated #7

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...

NO IMAGE
AVAILABLE

518134 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #5

6.2mm I.D., 9mm O.D., 68mm, fenestrated, inner and outer

NO IMAGE
AVAILABLE

506864 : FEDER-OSSOFF MICRO LARYNGEAL CUP FORCEPS

Angled up. Shaft tapers from 2.0 mm to 1 mm. Round biopsy cups w/ 1 mm jaw, 9" (23 cm)

NO IMAGE
AVAILABLE

507267S : ATOMIZER GLASS TUBE (SM #151)

NO IMAGE
AVAILABLE

507822 : JACKSON ESOPHAGEAL DILATOR

22 Fr.