

NO IMAGE
AVAILABLE

**518456 : JACKSON IMPROVED TRACHEOSTOMY TUBE #10,
W/15MM ADAPTER, REGULAR LENGTH**

10.6mm I.D., 14mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

507503 : JACKSON LARYNGEAL DILATOR

16 Fr., 11 3/4" (30cm)

NO IMAGE
AVAILABLE

511749 : TUCKER VALVE INNER TUBE TO FIT #8 TUCKER TUBE

9.1mm I.D., 63mm, sterling silver

NO IMAGE
AVAILABLE

**506860 : FEDER-OSSOFF MICRO LARYNGEAL GRASPING
FORCEPS**

Straight. Shaft tapers from 2.0 mm to 1 mm. Finely serrated, 9" (23 cm)

NO IMAGE
AVAILABLE

507505 : JACKSON LARYNGEAL DILATOR

20 Fr., 11 3/4" (30cm)

NO IMAGE
AVAILABLE

518115 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #2

3.4mm I.D., 6mm O.D., 51mm, stainless steel

NO IMAGE
AVAILABLE

**506858 : FEDER-OSSOFF MICRO LARYNGEAL GRASPING
FORCEPS**

Angled up. Shaft tapers from 2.0 mm to 1 mm. Finely serrated, 9" (23 cm)

NO IMAGE
AVAILABLE

507824 : JACKSON ESOPHAGEAL DILATOR

24 Fr.

NO IMAGE
AVAILABLE

518149 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #8

9.2mm I.D., 12mm O.D., 69mm, fenestrated, inner and outer

NO IMAGE
AVAILABLE

431325 : NOYES ALLIGATOR FORCEPS

Straight, Serrated, 8 x 1.8mm tip, 5 1/2" (14 cm)

NO IMAGE
AVAILABLE

507809 : JACKSON ESOPHAGEAL DILATOR

9 Fr.

NO IMAGE
AVAILABLE

**518335 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #7, SHORT
LENGTH**

8.1mm I.D., 11mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

507826 : JACKSON ESOPHAGEAL DILATOR

26 Fr.

NO IMAGE
AVAILABLE

**511125 : JACKSON IMPROVED TRACHEOSTOMY TUBE WITH
ROTATING LOCK #4, INNER AND OUTER TUBES AND
OBTURATOR**

5.3mm I.D., 8mm O.D., 62mm, sterling silver

NO IMAGE
AVAILABLE

431425 : Becker Scissors

Straight, Serrated, 7" (18.0cm)

NO IMAGE
AVAILABLE

505136 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Straight, 6mm cup, 40 cm

NO IMAGE
AVAILABLE

511744 : TUCKER VALVE INNER TUBE TO FIT #7 TUCKER TUBE

8.1mm I.D., 63mm, sterling silver

NO IMAGE
AVAILABLE

507500 : JACKSON LARYNGEAL DILATOR

10 Fr., 8 1/2" (21.5cm)

NO IMAGE
AVAILABLE

**518521 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #3, W/15MM
ADAPTER, REGULAR LENGTH**

4.3mm I.D., 7mm O.D., 58mm, stainless steel

NO IMAGE

**511140 : JACKSON IMPROVED TRACHEOSTOMY TUBE WITH
ROTATING LOCK #7, INNER AND OUTER TUBES AND
OBTURATOR**

8.1mm I.D., 11mm O.D., 69mm, sterling silver

AVAILABLE

NO IMAGE
AVAILABLE

507818 : JACKSON ESOPHAGEAL DILATOR

18 Fr.

NO IMAGE
AVAILABLE

511730 : TUCKER TUBE, INNER AND OUTER TUBES, SIZE 5

6.1mm I.D., 9mm O.D., 60mm, sterling silver

NO IMAGE
AVAILABLE

020204 : House Bellucci Alligator Scissors Mini

Straight, Sharp, 3 1/8" (8.0cm), 4.5mm blades

NO IMAGE
AVAILABLE

506465 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS

Straight, 8-5/8" (22 cm)

NO IMAGE
AVAILABLE

507504 : JACKSON LARYNGEAL DILATOR

18 Fr., 11 3/4" (30cm)