

NO IMAGE
AVAILABLE

507516 : JACKSON LARYNGEAL DILATOR

42 Fr., 11 3/4" (30cm)

NO IMAGE
AVAILABLE

073500 : Hurd Pillar Retractor, Small

7mm width x 1/4" (0.5cm) depth blade, 9" (23.0cm)

NO IMAGE
AVAILABLE

452551 : MAGILL CATHETER FORCEPS

Serrated, 7 1/8" (18 cm)

NO IMAGE
AVAILABLE

507515 : JACKSON LARYNGEAL DILATOR

40 Fr., 11 3/4" (30cm)

NO IMAGE
AVAILABLE

518129 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #4

5.5mm I.D., 8mm O.D., 62mm, fenestrated, inner and outer

NO IMAGE
AVAILABLE

1037826 : DECANNULATION STOPPER

Purple stopper with 40" cord, Size 6, 3/4 Closure

NO IMAGE
AVAILABLE

1037807 : DECANNULATION STOPPER

Yellow stopper with 40" cord, Size 7, 1/4 Closure

NO IMAGE
AVAILABLE

1037831 : DECANNULATION STOPPER

White stopper with 40" cord, Size 1, Full Closure

NO IMAGE
AVAILABLE

509203 : OSSOFF-KARLAN GRASPING FORCEPS EBONIZED

Straight, serrated, 2 mm tip, 8-3/8" (21 cm)

NO IMAGE
AVAILABLE

1037815 : DECANNULATION STOPPER

Green stopper with 40" cord, Size 5, 1/2 Closure

NO IMAGE
AVAILABLE

509006 : Angled Laser Retractor, Pediatric, Right

4mm diameter, 12 3/8" (31.5cm)

NO IMAGE
AVAILABLE

1037806 : DECANNULATION STOPPER

Yellow stopper with 40" cord, Size 6, 1/4 Closure

NO IMAGE
AVAILABLE

1037829 : DECANNULATION STOPPER

Purple stopper with 40" cord, Size 9, 3/4 Closure

NO IMAGE
AVAILABLE

505134 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Straight, 4mm cup, 60 cm

NO IMAGE
AVAILABLE

1037818 : DECANNULATION STOPPER

Green stopper with 40" cord, Size 8, 1/2 Closure

NO IMAGE
AVAILABLE

1037839 : DECANNULATION STOPPER

White stopper with 40" cord, Size 9, Full Closure

NO IMAGE
AVAILABLE

511144 : TUCKER VALVE INNER TUBE ONLY FOR #7 JACKSON IMPROVED TUBE

8.1mm I.D., 72mm, sterling silver

NO IMAGE
AVAILABLE

511130P : JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #5, INNER AND OUTER TUBES AND OBTURATOR

6.1mm I.D., 9.0mm O.D., 68mm, sterling silver

NO IMAGE
AVAILABLE

P26614 : PERKINS SELF-RETAINING RETRACTOR

3 blunt prongs on right, blade on left

NO IMAGE
AVAILABLE

505406 : JACKSON ALLIGATOR GRASPING FORCEPS

Straight, Serrated, 23 x 3.8mm tip, 23 5/8" (60 cm)

AVAILABLE

NO IMAGE
AVAILABLE

509010 : Laser Spatula Retractor, Pediatric, Right

5mm diameter, 22 cm working length

NO IMAGE
AVAILABLE

P26613 : PERKINS SELF-RETAINING RETRACTOR

3 blunt prongs on left, blade on right

NO IMAGE
AVAILABLE

507836 : JACKSON ESOPHAGEAL DILATOR

36 Fr.

NO IMAGE
AVAILABLE

525215 : Face Lift Retractor Fiberoptic Lighting On #2 Deaver Retractor

10-5/8" (27 cm)

NO IMAGE
AVAILABLE

20120 : ALLIGATOR FORCEPS

Serrated jaws, 6 mm, 3-1/8" (8 cm)