

NO IMAGE
AVAILABLE

514435 : Bernay Tracheal Retractor

Blunt, single prong, 9.5mm depth blade, 6 5/8" (17.0cm)

NO IMAGE
AVAILABLE

411525 : NOYES EAR FORCEPS

Serrated, 3 1/2" (9 cm)

NO IMAGE
AVAILABLE

031750 : Caplan Nasal Scissors

Angled 45, Serrated, 7" (18cm)

NO IMAGE
AVAILABLE

505140 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Straight, 6mm cup, 60 cm

NO IMAGE
AVAILABLE

070101 : Wieder Tongue Depressor

Fenestrated, serrated. Large, 44 mm blade. 6" (15 cm)

NO IMAGE
AVAILABLE

031800 : Knight Scissors

Angled, 42mm blades, 6 1/4" (16.0cm)

NO IMAGE
AVAILABLE

020210 : Bellucci Scissors

Curved right, Blunt, 3 3/8" (8.5cm)

NO IMAGE
AVAILABLE

506751 : JAKO FINE SPONGE CARRIER

22 cm

NO IMAGE
AVAILABLE

509008 : Laser Spatula Retractor, Adult, Right

10mm diameter, 22 cm working length

NO IMAGE
AVAILABLE

070116 : Jackson Trachea Retractor

Double-pronged, A-8.4mm / B-10.8mm width x A-8.5mm / B-16mm depth blade, 7" (18.0cm)

NO IMAGE
AVAILABLE

074140 : Dean Tonsil Scissors

Angled on flat, Serrated, 6 3/4" (17.0cm)

NO IMAGE
AVAILABLE

505146 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Angled, 4mm cup, 60 cm

NO IMAGE
AVAILABLE

506752 : JACKSON SPONGE CARRIER

Straight, 15-3/4" (40 cm)

NO IMAGE
AVAILABLE

505650 : PEDIATRIC BRONCHOSCOPIC GRASPING FORCEPS

Straight, Fine, Serrated, 14 1/8" (36 cm)

NO IMAGE
AVAILABLE

505674 : PILLING FLEXIBLE ENDOSCOPIC OVAL CUP FORCEPS

The Pilling Flexible Endoscopic Oval Cup Forceps were designed for use with Pilling Laser and Standard Bronchoscopes. They easily pass through the instrument/fiber channel for the purpose of biopsy, foreign body retrieval...

NO IMAGE
AVAILABLE

P27307 : Cottle Double Hook

Sharp prongs, 3mm wide, 5 7/8" (15cm)

NO IMAGE
AVAILABLE

509004 : Angled Laser Retractor, Adult, Right

8mm diameter, 12 3/4" (32.5cm)

NO IMAGE
AVAILABLE

514414 : Jackson Tracheal Retractor

Double-pronged, 3/8" (1.0cm) / 3/8" (1.0cm) width x 3/8"(1.0cm) / 1/2"(1.5cm) depth blade, 6 7/8" (17.5cm)

NO IMAGE
AVAILABLE

509201 : OSSOFF-KARLAN GRASPING FORCEPS EBONIZED

Curved left, 8-1/4" (21 cm)

NO IMAGE
AVAILABLE

505016 : ALLIGATOR GRASPING FORCEPS

Straight, Serrated, 2.4mm tip, 15" (38 cm)

AVAILABLE

NO IMAGE
AVAILABLE

505032 : KAHLER GRASPING FORCEPS

Straight, 19 5/8" (50 cm)

NO IMAGE
AVAILABLE

505148 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Angled, 6mm cup, 40 cm

NO IMAGE
AVAILABLE

509009 : Laser Spatula Retractor, Adult, Left

10mm diameter, 22 cm working length

NO IMAGE
AVAILABLE

P26646 : Bellucci Scissors

Straight 7mm blade, 5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

505113 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Straight, 2mm cup, 50 cm