

NO IMAGE
AVAILABLE

506075 : JACKSON FENESTRATED MEAT GRASPING FORCEPS

Straight, serrated, 19-5/8" (50 cm)

NO IMAGE
AVAILABLE

27100 : LUCAE BAYONET DRESSING FORCEPS

Serrated, 1.6 mm tip, 5-1/2" (14 cm)

NO IMAGE
AVAILABLE

411540 : HARTMAN-HERZFELD CUP FORCEPS

Straight, 3mm cups, 3 3/8" (9 cm)

NO IMAGE
AVAILABLE

514438 : Bernay Tracheal Retractor

Blunt, 4-pronged, 9.5mm depth blade, 6 5/8" (17.0cm)

NO IMAGE
AVAILABLE

P27663 : Foman Nasal Scissors

Curved, saber back blades, 5" (12.5cm)

NO IMAGE
AVAILABLE

P14311 : Knapp Retractor

6mm wide, 4 blunt prongs

NO IMAGE
AVAILABLE

20130 : HOUSE STRUT FORCEPS

Smooth jaws, 5 mm, 3-1/8" (8 cm)

NO IMAGE
AVAILABLE

P27598 : Sewell Orbital Retractor

No slot, 13mm x 67mm

NO IMAGE
AVAILABLE

P27599 : Sewell Orbital Retractor

No slot, 11mm x 50mm

NO IMAGE
AVAILABLE

505123 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Angled, 2mm cup, 50 cm

NO IMAGE
AVAILABLE

505152 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Angled, 6mm cup, 60 cm

NO IMAGE
AVAILABLE

411520 : HARTMANN DRESSING FORCEPS DELICATE

Fenestrated, Serrated, 2.5 x 7 mm jaws, 5 1/8" (13 cm)

NO IMAGE
AVAILABLE

P27681 : Lower Lateral Scissors

Fully curved blade, 6 1/2" (16.5cm)

NO IMAGE
AVAILABLE

P27630 : Converse Blade Retractor, #3

NO IMAGE
AVAILABLE

514437 : Bernay Tracheal Retractor

Blunt, 3-pronged, 9.5mm depth blade, 6 5/8" (17.0cm)

NO IMAGE
AVAILABLE

506750 : JACKSON SPONGE CARRIER

Straight, 11-7/8" (30 cm)

NO IMAGE
AVAILABLE

KM55170 : HOEN ALLIGATOR FORCEPS

Serrated, 1 x 2 teeth, 5-7/8" (15 cm)

NO IMAGE
AVAILABLE

P26647 : House Bellucci Alligator Scissors

Straight 5mm blade

NO IMAGE
AVAILABLE

31100 : LUCAE BAYONET DRESSING FORCEPS

Serrated, 2 mm tip, 6-1/2" (17 cm)

NO IMAGE

070107 : Andrews Depressor Large

20.5mm width blade

AVAILABLE

NO IMAGE
AVAILABLE

056360 : Fomon Lower Lateral Scissors

Strongly curved blades, 5 1/4" (13.3cm)

NO IMAGE
AVAILABLE

506076 : JACKSON FENESTRATED MEAT GRASPING FORCEPS

Straight, serrated, 23-5/8" (60 cm)

NO IMAGE
AVAILABLE

P26185 : Day Ear Hook, Large

90 deg.

NO IMAGE
AVAILABLE

P26183 : Day Ear Hook, Extra-Small

90 deg.

NO IMAGE
AVAILABLE

505048 : CIRCULAR CUP BIOPSY FORCEPS

Straight, 4mm cup, 45 cm