

NO IMAGE
AVAILABLE

KM55160 : HARTMAN ALLIGATOR FORCEPS

Serrated, 1.8 mm tip, 3-1/2" (9 cm)

NO IMAGE
AVAILABLE

505422 : JACKSON ALLIGATOR GRASPING FORCEPS

Straight, Serrated, 17 x 1.5mm tip, 15 3/4" (40 cm)

NO IMAGE
AVAILABLE

461605 : COLVER COAKLEY TONSIL FORCEPS

Medium Curved, 5 x 6 Teeth, 7 1/2" (19 cm)

NO IMAGE
AVAILABLE

505690 : PILLING FLEXIBLE ENDOSCOPIC CIRCULAR CUP FORCEPS

The Pilling Flexible Endoscopic Circular Cup Forceps were designed for use with Pilling Laser and Standard Bronchoscopes. They easily pass through the instrument/fiber channel for the purpose of biopsy, foreign body...

NO IMAGE
AVAILABLE

P27912 : BIOPSY AND GRASPING FORCEPS

3 mm cup, 70 deg., vertical

NO IMAGE
AVAILABLE

P27877 : WEIL BLAKESLEY FORCEPS, SIZE 0

Up 45 deg., 3.5 mm

NO IMAGE
AVAILABLE

505450 : JACKSON ROTATION GRASPING FORCEPS

Straight, Serrated, 4mm tip, 11 3/4" (30 cm)

NO IMAGE
AVAILABLE

505342 : JACKSON TRIANGULAR PUNCH FORCEPS

Angled, 3.5mm tip, 11 1/8" (28 cm)

NO IMAGE
AVAILABLE

P27851 : OSTRUM ANTRUM PUNCH FORCEPS

Right, 110 mm working length

NO IMAGE
AVAILABLE

P28177 : ST. CLAIRE THOMPSON ADENOID FORCEPS

7-1/4" (18.5 cm)

NO IMAGE
AVAILABLE

P28009 : THRUBITE BLAKESLEY FORCEPS, SIZE 0

Angled up, 3.5 mm

NO IMAGE
AVAILABLE

P27862 : OSTRUM ANTRUM PUNCH FORCEPS, PEDIATRIC

Rotating, 1.5 x 7 mm bite

NO IMAGE
AVAILABLE

P27882 : WEIL BLAKESLEY FORCEPS, SIZE 1

Straight, 4.0 mm



505138 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Straight, 6mm cup, 50 cm



20165 : ROUND CUP FORCEPS

0.75 mm cup jaws, 3-1/8" (8 cm)




P28001 : THRUBITE BLAKESLEY FORCEPS, SIZE 0

Straight, 3.5 mm




505150 : SAM ROBERTS LARYNGEAL CUP BIOPSY FORCEPS

Angled, 6mm cup, 50 cm




452510 : COLLINGS TONGUE SEIZING FORCEPS

Straight, Serrated, 6 3/4" (17 cm)



443610 : ASCH SEPTUM FORCEPS

Angled, 28mm tip, 9 1/8" (23 cm)



57150 : ALLIGATOR NASAL FORCEPS

Serrated jaw, 10 x 1 mm, 7" (18 cm)

AVAILABLE

NO IMAGE
AVAILABLE

506070 : JACKSON PEANUT GRASPING FORCEPS

Straight, serrated, 23-5/8" (60 cm)

NO IMAGE
AVAILABLE

505680 : PILLING FLEXIBLE ENDOSCOPIC FORCEPS

Straight, Serrated, 1.5mm tip, 15 3/4" (40 cm)

NO IMAGE
AVAILABLE

P27878 : WEIL BLAKESLEY FORCEPS, SIZE 1

Up 45 deg., 4.0 mm

NO IMAGE
AVAILABLE

P27880 : WEIL BLAKESLEY FORCEPS, SIZE 0

Straight, 3.5 mm

NO IMAGE
AVAILABLE

505160 : BIOPSY CUP FORCEPS

Straight, 2mm cup, 13 3/4" (35 cm)