**461600 : COLVER COAKLEY TONSIL FORCEPS**

Straight, 5 x 6 Teeth, 7 1/2" (19 cm)

**461660 : MELTZER ADENOID PUNCH SIZE 1**

Straight, 5mm, Working length 5 1/2" (14cm), Size 1

**461661 : MELTZER ADENOID PUNCH SIZE 2**

Straight, 7mm, Working length 5 1/2" (14cm), Size 2

**502245 : Laryngoscope Holder**

Pilling Laryngoscope Holders enable the surgeon to have both hands free during Laryngoscopy, Micro-Laryngoscopy and Laser Laryngoscopy procedures. These holders allow the surgeon to perform intricate procedures requiring...

**502248 : Tobey Laryngoscope Holder**

Pilling Laryngoscope Holders enable the surgeon to have both hands free during Laryngoscopy, Micro-Laryngoscopy and Laser Laryngoscopy procedures. These holders allow the surgeon to perform intricate procedures requiring...

**502260 : Boston University Suspension System**

This system is used to suspend the laryngoscope during micro and laser surgery of the larynx. *Laryngoscope and Strong Laryngoscope Handle Extension (502261) are not included.* • Supports the laryngoscope for two-handed...



### **502261 : Laryngoscope Handle Extension for Boston University System**

Adapts micro-laryngoscope for use with the Boston University Suspension System. Laryngoscopes already in use may be returned to Teleflex Instrument Repair for the modification or a new instrument may be ordered with this...



### **502262 : Boston University System Pivoting Extension Arm Assembly**

502262: BU Suspension Extension Arm, overall length 10" (25.5 cm)



### **502275 : Teeth Guards**

These nylon teeth guards were designed as protection for the patient's teeth against pressure that may be exerted during insertion and positioning of a laryngoscope. 502275: Teeth Guards Adult size teeth guards, box of 24...



### **502276 : Teeth Guards**

These nylon teeth guards were designed as protection for the patient's teeth against pressure that may be exerted during insertion and positioning of a laryngoscope. 502276: Teeth Guards Adult size teeth guards, box of 3...



### **502540 : FLUVOG SLIDE ASSEMBLY FOR PILLING BRONCHOSCOPES**



### **502566 : Sanders Venturi Adapter For Pilling® Bronchoscopes**

Flow and pressure characteristics vary greatly according to the design of the bronchoscope or laryngoscope and any new combination should be studied carefully before clinical use. *Note: The Sanders Venturi System can be...*



### **502594 : VENTURI VENTILATING NEEDLE**

14 gauge tip, 5-3/8" (13.8 cm)



**502599 : Replacement Luer-Lok® Hub Adapter**

Needle Luer-Lok to Sanders Tubing, 1-1/8" (3 cm)



**502600 : VENTURI VENTILATING NEEDLE**

Replacement needle for 522244, 522245, 522242, and 522229. 14 gauge tip, 4" (10.2cm)



**505010 : FENESTRATED CUPPED GRASPING FORCEPS**

Straight, 3.8mm tip, 15" (38 cm)



**505022 : ROTATION GRASPING FORCEPS**

Straight, Serrated, 1 x 1 Teeth, 15" (38 cm)



**505050 : CIRCULAR CUP BIOPSY FORCEPS**

Straight, 4mm cup, 50 cm




**505334 : HOLINGER CUP BITING SPECIMEN FORCEPS**

Straight, 23 5/8" (60 cm)




**505455 : JACKSON ROTATION GRASPING FORCEPS**

Straight, Serrated, 4mm tip, 19 5/8" (50 cm)




**505456 : JACKSON ROTATION GRASPING FORCEPS**

Straight, Serrated, 4mm tip, 23 5/8" (60 cm)




**505682 : PILLING FLEXIBLE ENDOSCOPIC FORCEPS**

Straight, Serrated, 2.2mm tip, 15 3/4" (40 cm)




**505686 : PILLING FLEXIBLE ENDOSCOPIC FORCEPS**

Straight, Serrated, 19 5/8" (50 cm)



**506003 : JACKSON FORWARD GRASPING CUP FORCEPS**

Straight, serrated, 13-3/4" (35 cm)



**506004 : JACKSON FORWARD GRASPING CUP FORCEPS**

Straight, serrated, 15-3/4" (40 cm)