

Home // ENT

ENT

**506830 : FEDER-OSSOFF PROBE**

30 angle, blunt, 10 7/8" (27.5cm)

**506832 : FEDER-OSSOFF PROBE**


60 angle, blunt, 10 5/8" (27cm)

**506842 : FEDER-OSSOFF PROBE**


30 angle, sharp, 10 5/8" (27cm)

**506844 : FEDER-OSSOFF PROBE**


30 angle, sharp, 10 5/8" (27cm)

**506870 : FEDER-OSSOFF SUCTION TUBE**

3 Fr. tip

**506872 : FEDER-OSSOFF SUCTION TUBE**

5 Fr. tip



506874 : FEDER-OSSOFF SUCTION TUBE

7 Fr. tip



506880 : FEDER-OSSOFF SICKLE KNIFE

Straight sickle, 3mm blade, 27cm working length, 3.2mm - 1mm tapered shaft.



506882 : FEDER-OSSOFF SICKLE KNIFE

Curved sickle, 3mm blade, 27cm working length, 3.2mm - 1mm tapered shaft.



506884 : FEDER-OSSOFF SICKLE KNIFE

Phono circle knife, 3mm blade, 27cm working length, 3.2mm - 1mm tapered shaft.



506888 : Feder-Ossoff Universal Ring Handled Knife Handle

With camlock506888: Feder-Ossoff Universal Ring Handled Knife Handle



506890 : FEDER-OSSOFF PHONOSURGERY STERILIZING TRAY



507010 : JACKSON OPEN END ASPIRATING TUBE

12 Fr. tip, working length 7-7/8" (20 cm)

507011 : JACKSON OPEN END ASPIRATING TUBE

12 Fr. tip, working length 13-3/4" (35 cm)



507012 : JACKSON OPEN END ASPIRATING TUBE

12 Fr. tip, working length 15-3/4" (40 cm)



507013 : JACKSON OPEN END ASPIRATING TUBE

12 Fr. tip, working length 17-3/4" (45 cm)



507014 : JACKSON OPEN END ASPIRATING TUBE

12 Fr. tip, working length 19-5/8" (50 cm)



507015 : JACKSON OPEN END ASPIRATING TUBE

12 Fr. tip, working length 23-5/8" (60 cm)



507020 : JACKSON OPEN END ASPIRATING TUBE


6 Fr. tip, working length 13-3/4" (35 cm)



507023 : JACKSON OPEN END ASPIRATING TUBE

8 Fr. tip, working length 13-3/4" (35 cm)





507027 : JACKSON OPEN END ASPIRATING TUBE

9 Fr. tip, working length 13-3/4" (35 cm)



507028 : JACKSON OPEN END ASPIRATING TUBE

9 Fr. tip, working length 15-3/4" (40 cm)



507030 : JACKSON VELVET EYE ASPIRATING TUBE

12 Fr. tip, working length 7-7/8" (20 cm)



507031 : JACKSON VELVET EYE ASPIRATING TUBE

12 Fr. tip, working length 13-3/4" (35 cm)



507032 : JACKSON VELVET EYE ASPIRATING TUBE

12 Fr. tip, working length 15-3/4" (40 cm)

