**507033 : JACKSON VELVET EYE ASPIRATING TUBE**

12 Fr. tip, working length 17-3/4" (45 cm)

**507034 : JACKSON VELVET EYE ASPIRATING TUBE**


12 Fr. tip, working length 19-5/8" (50 cm)

**507035 : JACKSON VELVET EYE ASPIRATING TUBE**

12 Fr. tip, working length 23-5/8" (60 cm)

**507040 : JACKSON VELVET EYE ASPIRATING TUBE**

9 Fr. tip, working length 13-3/4" (35 cm)

**507041 : JACKSON VELVET EYE ASPIRATING TUBE**

12 Fr. tip, working length 15-3/4" (40 cm)

**507045 : INSULATED SUCTION CAUTERY TUBE**

12 Fr. tip, working length 8-5/8" (22 cm), supplied w/10 ft. Bovie cord and plug (507046)



507046 : BOVIE CORD FOR SUCTION TUBE

Replacement bovie cord and plug for INSULATED SUCTION CAUTERY TUBE (507045) and FULLY INSULATED SUCTION TUBE (507048), 10 ft.



507048 : FULLY INSULATED SUCTION TUBE

12 Fr. tip, working length 19-5/8" (50 cm), supplied w/10 ft. Bovie cord and plug (507046)



507050 : HOLINGER ASPIRATING TUBE

12 Fr. tip, working length 17-3/4" (45 cm)



507051 : HOLINGER ASPIRATING TUBE

12 Fr. tip, working length 19-5/8" (50 cm)



507052 : HOLINGER ASPIRATING TUBE

12 Fr. tip, working length 23-5/8" (60 cm)



507055 : HOLINGER ASPIRATING TUBE

9 Fr. tip, working length 9-7/8" (25 cm)



507057 : HOLINGER ASPIRATING TUBE

9 Fr. tip, working length 13-3/4" (35 cm)




507060 : HOLINGER ASPIRATING TUBE

9 Fr. tip, working length 7-7/8" (20 cm)



507065 : HOLINGER ASPIRATING TUBE

12 Fr. tip, working length 13-3/4" (35 cm)




507120 : JAKO SUCTION TUBE SIZE 00

5 Fr. tip




507121 : JAKO SUCTION TUBE SIZE 0

6 Fr. tip



507122 : JAKO SUCTION TUBE SIZE 1

7 Fr. tip



507453 : ASPIRATING NEEDLE

Reinforced angular shaft, with protective cover, distal needle 23 gauge (0.6 mm), 14.7 mm long, working length 8 1/4" (21 cm)



507454 : SPARE DISTAL NEEDLE TIP FOR ASPIRATING NEEDLE 507453

23 gauge (0.6 mm), 14.7 mm long


507498S : ATOMIZER REPLACEMENT BULB (SM #1-RD)

6 1/8" (15.5cm)



507807 : JACKSON ESOPHAGEAL DILATOR SET

Includes nineteen dilators, 8 Fr. - 42 Fr.



508200 : LARYNGEAL MIRROR

#000, 8mm diameter, Shaft length is 6 3/4" (17.0cm)



508201 : LARYNGEAL MIRROR

#00, 10mm diameter, Shaft length is 6 3/4" (17.0cm)



508202 : LARYNGEAL MIRROR

#0, 12mm diameter, Shaft length is 6 3/4" (17.0cm)

