**508203 : LARYNGEAL MIRROR**


#1, 16mm diameter, Shaft length is 6 3/4" (17.0cm)

**508204 : LARYNGEAL MIRROR**

#2, 18mm diameter, Shaft length is 6 3/4" (17.0cm)

**508205 : LARYNGEAL MIRROR**

#3, 20mm diameter, Shaft length is 6 3/4" (17.0cm)

**508206 : LARYNGEAL MIRROR**


#4, 22mm diameter, Shaft length is 6 3/4" (17.0cm)

**508207 : LARYNGEAL MIRROR**

#5, 24mm diameter, Shaft length is 6 3/4" (17.0cm)

**508208 : LARYNGEAL MIRROR**

#6, 28mm diameter, Shaft length is 6 3/4" (17.0cm)



508209 : LARYNGEAL MIRROR

#7, 30mm diameter, Shaft length is 6 3/4" (17.0cm)




508230 : JAKO MIRROR

Working length 8 1/4" (21.0cm)



508540 : Holman Brushes

3/set




508545 : Trachea Tube Brushes

3/Set



508590 : JACKSON FORCEPS PARTS KIT



509012 : JAKO VELVET EYE ASPIRATING TUBE

Working length 7 7/8" (20cm)



509013 : JAKO ASPIRATING TUBE-OPEN END

Open end, working length 6" (15cm)

509014 : ASPIRATING TUBE CLIP

Clip for attaching aspirating tubes to laryngoscope up to 9 Fr. (3 mm)




509204 : OSSOFF MICRO SUCTION TUBE EBONIZED

7 7/8" (20.0cm)



509205 : MICRO-LARYNGEAL PROBE

Straight, Ebonized, 10 1/4" (26cm)



509206 : MICRO-LARYNGEAL PROBE

Angled, Ebonized, 10" (25.5cm)




509209 : LASER MIRROR

Highly polished stainless front surface, 6mm diameter, Working length 8 1/4" (21.0cm)



509301 : MICROSUBGLOTIC OPERATING PLATFORM 3MM X 28CM

0.1" (3mm), 14 1/2" (37.0cm)



509302 : FRONT SURFACE MIRROR

3mm diameter, Working length 11" (28.0cm)



509304 : MICROSUBGLOTTIC STRAIGHT PROBE

Straight, 13" (33cm)



509305 : SUBGLOTTIC GRASPING FORCEPS

Straight, serrated, 11" (28 cm)



511149 : TUCKER VALVE INNER TUBE ONLY FOR #8 JACKSON IMPROVED TUBE

9.1mm I.D., 72mm



511745 : TUCKER TUBE, INNER AND OUTER TUBES, OBTURATOR, SIZE 8

9.1mm I.D., 12mm O.D., 60mm, sterling silver



511930 : MAYO CLINIC TUBE #5

6.2mm I.D., 9mm O.D., 56mm, stainless steel

