

**511940 : MAYO CLINIC TUBE #6**

7.2mm I.D., 10mm O.D., 56mm, stainless steel

**511945 : MAYO CLINIC TUBE #6**

7.2mm I.D., 10mm O.D., 69mm, stainless steel

**511955 : MAYO CLINIC TUBE #7**

8.1mm I.D., 11mm O.D., 69mm, stainless steel

**512780 : NORRIS TRACHEOSTOMY BUTTONS**

8.3mm I.D., 19mm x 11mm, white teflon

**514436 : Bernay Tracheal Retractor**

Blunt, Double-pronged, 9.5mm depth blade, 6 5/8" (17.0cm)

**516270 : HOLINGER LARYNGECTOMY TUBE #7**

9.0mm I.D., 12mm O.D., 50mm, sterling silver



**516275 : HOLINGER LARYNGECTOMY TUBE #8**

10.0mm I.D., 13.6mm O.D., 50mm, sterling silver



**516280 : HOLINGER LARYNGECTOMY TUBE #9**

10.6mm I.D., 14mm O.D., 50mm, sterling silver



**516285 : HOLINGER LARYNGECTOMY TUBE #10**

11.2mm I.D., 14.3mm O.D., 50mm, sterling silver



**516430 : MARTIN LARYNGECTOMY TUBE, EXTENSION SET, SIZE 12**

12.7mm I.D., 16.3mm O.D., 42mm, sterling silver




**518025A : JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #4 FOR SHIKANI-FRENCH SPEAKING VALVE**

5.5mm I.D., 8mm O.D., 62mm



**518027 : Non-Fenestrated Jackson Improved Inner Tube Only with 2 Shikani Valves, Size 4**

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...



**518029A : JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #4 FOR SHIKANI-FRENCH SPEAKING VALVE**

5.5mm I.D., 8mm O.D., 62mm, fenestrated, inner and outer

NO IMAGE  
AVAILABLE

**518034A : JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #5 FOR SHIKANI-FRENCH SPEAKING VALVE**

6.2mm I.D., 9.0mm O.D., 68mm, fenestrated, inner and outer

**518037 : Non-Fenestrated Jackson Improved Inner Tube Only with 2 Shikani Valves, Size 6**

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...

**518042 : Non-Fenestrated Jackson Improved Inner Tube Only with 2 Shikani Valves, Size 7**

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...

**518315 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #3, SHORT LENGTH**

4.3mm I.D., 7mm O.D., 48mm, stainless steel

**518556 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #10, W/15MM ADAPTER, REGULAR LENGTH**

10.6mm I.D., 14mm O.D., 69mm, stainless steel

**518640 : JACKSON IMPROVED TRACHEOSTOMY TUBE #8, W/15MM ADAPTER, SHORT LENGTH**

9.2mm I.D., 12mm O.D., 56mm, stainless steel

**521450 : TUCKER MID-LIGHTED SIDE SLIDE LARYNGOSCOPE**

Overall length 91mm, min. I.D. 6.9mm, replacement light carrier 526510



**521455 : TUCKER MID-LIGHTED SIDE SLIDE LARYNGOSCOPE, ADOLESCENT**

Overall length 106mm, min. I.D. 9.0mm, replacement light carrier 526512



**521460 : TUCKER MID-LIGHTED SIDE SLIDE LARYNGOSCOPE, ADULT**

Overall length 121mm, min. I.D. 10.7mm, replacement light carrier 526515



**521465 : TUCKER MID-LIGHTED SIDE SLIDE LARYNGOSCOPE, ADULT LONG**

Overall length 160mm, min. I.D. 11.8mm, replacement light carrier 526516



**521500 : HOLINGER LARYNGOSCOPE, INFANT**

Overall length 82mm, min. I.D. 8.0mm, replacement light carrier 526204



**521550 : JACKSON SLIDE LARYNGOSCOPE, INFANT**

Overall length 95mm, min. I.D. 8.0mm, replacement light carrier 526000

