

**521552 : JACKSON SLIDE LARYNGOSCOPE, CHILD**

Overall length 125mm, min. I.D. 10.5mm, replacement light carrier 526006

**521554 : JACKSON SLIDE LARYNGOSCOPE, ADOLESCENT**

Overall length 160mm, min. I.D. 10.6mm, replacement light carrier 526018

**521556 : JACKSON SLIDE LARYNGOSCOPE, ADULT**

Overall length 160mm, min. I.D. 11.4mm, replacement light carrier 526018

**521558 : JACKSON SLIDE LARYNGOSCOPE, ADULT**

Overall length 160mm, min. I.D. 13.0mm, replacement light carrier 526018

**521820 : JACKSON ANTERIOR COMMISSURE LARYNGOSCOPE, CHILD**

Overall length 133mm, min. I.D. 5.8mm, replacement light carrier 526008

**521824 : JACKSON ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT**

Overall length 173mm, min. I.D. 6.8mm, replacement light carrier 526026



522002 : HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, INFANT

Overall length 97mm, replacement light carrier 526206



522004 : HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, CHILD

Overall length 120mm, replacement light carrier 526210



522010 : HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT

Overall length 170mm, min. I.D. 8.3mm, replacement light carrier 526028



522012 : HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT EXTRA-LONG

Overall length 180mm, min. I.D. 8.3mm, replacement light carrier 526038



522030 : HOLINGER "HOUR GLASS" ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT

Overall length 171mm, min. I.D. 7.7mm, replacement light carrier 526212



522032 : HOLINGER "HOUR GLASS" ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT EXTRA LONG

Overall length 181mm, min. I.D. 7.3mm, replacement light carrier 526214



522040 : HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT SLOTTED

Overall length 170mm, min. I.D. 8.5mm, replacement light carrier 526028



522045 : HOLINGER "HOUR GLASS" ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT SLOTTED

Overall length 170mm, min. I.D. 8.2mm, replacement light carrier 526212



522190 : Ossoff-Pilling® Micro-Laryngoscope, Adult Female

The new Ossoff-Pilling Micro-Laryngoscope is available in two proximal widths. The 25.5 mm width (522191) can be introduced in most patients and is easily aligned with the microscope (400 mm objective). The 19 mm scope...



522191 : Ossoff-Pilling® Micro-Laryngoscope, Adult Male

The new Ossoff-Pilling Micro-Laryngoscope is available in two proximal widths. The 25.5 mm width (522191) can be introduced in most patients and is easily aligned with the microscope (400 mm objective). The 19 mm scope...




522195 : SUCTION TUBE FOR 522194 AND 522226

7 1/8" (18.0cm)



522196 : Garrett-Ossoff-Pilling® Laryngoscope, Adult Extra-long

The new Garrett-Ossoff-Pilling Laryngoscope helps surgeons focus on critical anatomy during procedures in a range of patient sizes, and is designed for patients with a greater distance between their teeth and vocal cords....



522200P : JAKO LARYNGOSCOPE TYPE 1 ADULT

NO IMAGE AVAILABLE



522205 : JAKO LARYNGOSCOPE TYPE 2 ADULT

NO IMAGE

AVAILABLE



522206 : JAKO MICRO LARYNGOSCOPE TYPE III, ADULT

Overall length 178mm, min. I.D. 13.4mm, replacement light carrier 526526



522207 : JAKO-CHERRY LARYN TYPE 1 ADOL



522208 : JAKO-CHERRY LARYN TYPE 2 ADOL



522209 : JAKO MICRO LARYNGOSCOPE TYPE III, ADOLESCENT

Overall length 149mm, min. I.D. 11.7mm, replacement light carrier 526527



522216 : JAKO MICRO LARYNGOSCOPE TYPE III, CHILD

Overall length 120mm, min. I.D. 9.1mm, replacement light carrier 526528