

**522220 : JAKO-CHERRY MICRO LARYNGOSCOPE TYPE II, NEWBORN**

Small, Overall length 114mm, min. I.D. 7.8mm, replacement light carrier 527825, 527820

**522221 : Low-Profile Dedo-Style Laser Laryngoscope, Adult**

The Newest Addition to the Pilling Family. Featuring a low-profile design to accommodate a wider range of patient anatomies than previous-generation Pilling Dedo laryngoscopes, while maintaining the same internal diameter...

**522222 : Dedo-Pilling® Micro-Laryngoscope, Child**


Our Most Popular. The Dedo-Pilling Micro-Laryngoscopes have been designed to provide the ability to examine the vocal cords at the anterior commissure in the majority of patients, for example, patients with bull necks,...

**522224 : Dedo-Pilling® Micro-Laryngoscope, Adolescent**

Our Most Popular. The Dedo-Pilling Micro-Laryngoscopes have been designed to provide the ability to examine the vocal cords at the anterior commissure in the majority of patients, for example, patients with bull necks,...

**522225P : Dedo-Pilling® Micro-Laryngoscope, Adult**

Our Most Popular. The Dedo-Pilling Micro-Laryngoscopes have been designed to provide the ability to examine the vocal cords at the anterior commissure in the majority of patients, for example, patients with bull necks,...

**522227 : DEDO MICRO LARYNGOSCOPE TYPE 2**




**522228 : DEDO LASER MICRO LARYNGOSCOPE, ADULT**

Overall length 177mm, min. I.D. 10.7mm, replacement light carrier 527860



**522229 : ABRAMSON-DEDO MICRO LARYNGOSCOPE, ADULT**

Overall length 177mm, min. I.D. 12.0mm, replacement light carrier 527586, 527587



**522231 : OSSOFF-KARLAN-DEDO LASER LARYNGOSCOPE, ADULT**

Overall length 178mm, min. I.D. 11.9mm, replacement light carrier 527860



**522232 : OSSOFF-KARLAN-JAKO LASER LARYNGOSCOPE**

Overall length 178mm, min. I.D. 13.4mm, replacement light carrier 527800, 527805



**522233 : OSSOFF-KARLAN POSTERIOR COMMISSURE LARYNGOSCOPE, ADULT**

Overall length 178mm, min. I.D. 11.9mm, replacement light carrier 527860



**522237 : OSSOFF-KARLAN LASER LARYNGOSCOPE, CHILD**

Overall length 130mm, min. I.D. 9.8mm, replacement light carrier 527863



**522238 : OSSOFF-KARLAN LASER LARYNGOSCOPE, ADOLESCENT**

Overall length 152mm, min. I.D. 9.7mm, replacement light carrier 527862



**522240 : REPLACEMENT CLIP FOR VENTURI VENTILATING NEEDLE ON 522229**

5/8" (1.6cm)



**522241 : VENTURI VENTILATING NEEDLE**

Blunt, 14 gauge, bent 90, 12" (30.5cm)



**522242 : HEALY-JAKO SUBGLOTTISCOPE, PEDIATRIC**

Overall length 155mm, min. I.D. 5.7mm, replacement light carrier 527830 right, 527835 left



**522243 : DEDO-PILLING LASER LARYNGOSCOPE, ADULT**

Large, Overall length 197mm, min. I.D. 12.6mm, replacement light carrier 527860



**522244 : OSSOFF SUBGLOTTISCOPE, FEMALE ADULT**

Large, Overall length 230mm, min. I.D. 7.3mm, replacement light carrier 527860



**522245 : OSSOFF SUBGLOTTISCOPE, MALE ADULT**

Large, Overall length 230mm, min. I.D. 9.5mm, replacement light carrier 527860



**522290 : PILLING-DEDO DOUBLE-ACTION ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT**

Large, min. I.D. 12.7mm, replacement light carrier 527860



**522291 : PILLING-OSSOFF DOUBLE-ACTION POSTERIOR COMMISSURE LARYNGOSCOPE, ADULT**

Large, min. I.D. 12.7mm, replacement light carrier 527860



**522292 : PILLING-DEDO DOUBLE-ACTION ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT**

min. I.D. 9.8mm, replacement light carrier 527863



**522550 : PILLING-HOLINGER BRONCHOSCOPE 3X20**

Infant, Scope size 3 x 20, overall length 203mm, min. I.D. 2.9mm, replacement light carrier 527500



**522552 : PILLING-HOLINGER BRONCHOSCOPE 3.5X20**

Infant, Scope size 3.5 x 20, overall length 203mm, min. I.D. 3.5mm, replacement light carrier 527500



**522554 : PILLING-HOLINGER BRONCHOSCOPE 4X20**

Infant, Scope size 4 x 20, overall length 203mm, min. I.D. 3.9mm, replacement light carrier 527500

