

**522556 : PILLING-HOLINGER BRONCHOSCOPE 4X25**

Infant, Scope size 4 x 25, overall length 254mm, min. I.D. 3.9mm, replacement light carrier 527505

**522575 : HOLINGER BRONCHOSCOPE 3X20**


Scope size 3 x 20, overall length 203mm, min. I.D. 2.5mm, replacement light carrier 527500

**522576 : HOLINGER BRONCHOSCOPE 3.5X25**


Scope size 3.5 x 25, overall length 254mm, min. I.D. 3.4mm, replacement light carrier 527505

**522577 : HOLINGER BRONCHOSCOPE 3.5X30**

Scope size 3.5 x 30, overall length 305mm, min. I.D. 3.4mm, replacement light carrier 527510

**522578 : HOLINGER BRONCHOSCOPE 4X30**

Scope size 4 x 30, overall length 305mm, min. I.D. 3.7mm, replacement light carrier 527510

**522579 : HOLINGER BRONCHOSCOPE 5X30**

Scope size 5 x 30, overall length 305mm, min. I.D. 4.5mm, replacement light carrier 527510



**522580 : HOLINGER BRONCHOSCOPE 5X35**

Scope size 5 x 35, overall length 350mm, min. I.D. 4.8mm, replacement light carrier 527515



**522581 : HOLINGER BRONCHOSCOPE 5X40**

Scope size 5 x 40, overall length 400mm, min. I.D. 4.8mm, replacement light carrier 527520



**522582 : HOLINGER BRONCHOSCOPE 5X45**

Scope size 5 x 45, overall length 451mm, min. I.D. 4.9mm, replacement light carrier 527525



**522583 : HOLINGER BRONCHOSCOPE 6X35**

Scope size 6 x 35, overall length 350mm, min. I.D. 5.6mm, replacement light carrier 527555



**522584 : HOLINGER BRONCHOSCOPE 6X40**

Scope size 6 x 40, overall length 400mm, min. I.D. 5.6mm, replacement light carrier 527560



**522585 : HOLINGER BRONCHOSCOPE 7X40**

Scope size 7 x 40, overall length 400mm, min. I.D. 5.8mm, replacement light carrier 527560



**522586 : HOLINGER BRONCHOSCOPE 8X40**

Scope size 8 x 40, overall length 400mm, min. I.D. 7.6mm, replacement light carrier 527580

**522587 : HOLINGER BRONCHOSCOPE 9X40**

Scope size 9 x 40, overall length 400mm, min. I.D. 8.5mm, replacement light carrier 527580



**523120P : LASER BRONCHOSCOPE 5X30**

Scope size 5.0 x 30, overall length 300mm, min. I.D. 5.1mm, replacement light carrier 527550



**523130P : LASER BRONCHOSCOPE 6X35**

Scope size 6.0 x 35, overall length 350mm, min. I.D. 6.0mm, replacement light carrier 523131



**523131 : LIGHT CARRIER OVAL LARGE 35CM**

For scope number 523130P



**523140P : LASER BRONCHOSCOPE 7X40**

Scope size 7.0 x 40, overall length 400mm, min. I.D. 7.0mm, replacement light carrier 527580



**523150P : LASER BRONCHOSCOPE 8X40**


Scope size 8.0 x 40, overall length 400mm, min. I.D. 8.0mm, replacement light carrier 527580



**523160P : LASER BRONCHOSCOPE 9X40**

Scope size 9.0 x 40, overall length 400mm, min. I.D. 9.0mm, replacement light carrier 527580





**523210P : MOISTURE BARRIER ASSEMBLY FOR CO2 COUPLER**

1 5/8" (4.0cm)



**523505 : HOLINGER ESOPHAGOSCOPE**

Infant, 6 x 20, Overall length 200mm, min. I.D. 4.0mm, replacement light carrier 526218



**523570 : HOLINGER ESOPHAGOSCOPE**

Infant, 3.5 x 30, Overall length 300mm, min. I.D. 3.4mm, replacement light carrier 527510



**523571 : HOLINGER ESOPHAGOSCOPE**

4 x 30, Overall length 300mm, min. I.D. 3.9mm, replacement light carrier 527510



**523572 : HOLINGER ESOPHAGOSCOPE**

6 x 30, Overall length 300mm, min. I.D. 5.4mm, replacement light carrier 527550