

NO IMAGE  
AVAILABLE

**528308P : XenaLight™ XLT Bulb with Timer Board for Older Model  
528305, 528305XLT**

NO IMAGE  
AVAILABLE

**528309P : XenaLight™ XLT II Lighting System FLOORSTAND  
Base/Pole/Swivel**

NO IMAGE  
AVAILABLE

**528310P : XenaLight™ XLT II Bulb with Timer Chip for New Model  
528305XLTP**

NO IMAGE  
AVAILABLE

**528339P : Variable 10-75 mm Spot SURGILUX™ HEADLIGHT  
MODULE with 8' XLG BIF Cable**

NO IMAGE  
AVAILABLE

**528340P : Fixed 100 mm Spot SURGILUX™ HEADLIGHT MODULE  
with 8' XLG BIF Cable**

NO IMAGE  
AVAILABLE

**528341P : Variable 10-75 mm Spot SURGILUX™ HEADLIGHT  
MODULE with 8' XLG Cable**

NO IMAGE  
AVAILABLE

**528342P : Fixed 100 mm Spot SURGILUX™ HEADLIGHT MODULE  
with 8' XLG Cable**

NO IMAGE  
AVAILABLE

**528343P : Variable 10-75 mm Spot SURGILUX™ HEADLIGHT  
MODULE with 10' XLG Cable**

NO IMAGE  
AVAILABLE

**528344P : Fixed 100 mm Spot SURGILUX™ HEADLIGHT MODULE  
with 10' XLG Cable**

NO IMAGE  
AVAILABLE

**528351P : Headband Front Insert, Breath-O-Prene (Pkg of 5)**

NO IMAGE  
AVAILABLE

**528352P : Headband Top Insert, Breath-O-Prene (Pkg of 5)**

NO IMAGE  
AVAILABLE

**528353P : Headband Back Insert, Breath-O-Prene (Pkg of 5)**

NO IMAGE  
AVAILABLE

**528354P : Joysticks (Pkg of 5)**



**528355P : Gown Clips (Pkg of 10)**



**70501 : La Force Adenotome Size 1**

9" (23 cm). No. 1, small, blade 25 x 16 mm




**70502 : La Force Adenotome Size 2**

9" (23 cm). No. 2, medium, blade 26 x 19 mm




**70503 : La Force Adenotome Size 3**

9" (23 cm). No. 3, large, blade 30 x 22 mm




**70511 : Shambaugh Reverse Curve Adenotomes**

8-5/8" (22 cm). No. 1, small, blade 25 x 16 mm




**70512 : Shambaugh Reverse Curve Adenotomes**

8-5/8" (22 cm). No. 2, medium, blade 26 x 19 mm



**A358AFA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. A358AFA: Pilling Blue Fiber Optic Cable, ACMI illuminator end, 3.5 mm bundle, 8 ft. length, ACMI Female instrument end



AVAILABLE

NO IMAGE  
AVAILABLE

**A358PRA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. A358PRA: Pilling Blue Fiber Optic Cable, ACMI illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end (formally A358PR)

NO IMAGE  
AVAILABLE

**A358STA : Fiberoptic Cable**

ACMI illuminator end, 3.5 mm bundle, 8' length, Storz threaded instrument end

NO IMAGE  
AVAILABLE

**A508STA : Fiberoptic Cable**

ACMI illuminator end, 5 mm bundle, 8' length, Storz threaded instrument end

NO IMAGE  
AVAILABLE

**KM55094 : LUCAE BAYONET DRESSING FORCEPS**

Serrated, 1.8 mm tip, 5-1/2" (14 cm)

NO IMAGE  
AVAILABLE

**KM55100 : JANSEN BAYONET DRESSONG FORCEPS**

Serrated, 2.1 mm tip, 6-1/4" (16 cm)