

NO IMAGE
AVAILABLE

P26395 : LANCET KNIFE

Angled 45 deg., 2 x 3mm tip, 6 1/2" (16.5cm)

NO IMAGE
AVAILABLE

P26396 : SHEEHY KNIFE CURETTE

Semi-sharp, 2mm dia. blade

NO IMAGE
AVAILABLE

P26397 : SHEEHY KNIFE CURETTE

Semi-sharp, 2.5mm dia. blade

NO IMAGE
AVAILABLE

P26398 : TYMPANOPLASTY KNIFE

Curved blade, 7mm, 6 1/2" (16.5cm)

NO IMAGE
AVAILABLE

P26399 : INCUDOSTAPEDIAL JOINT KNIFE

Angled blade, 0.75mm dia.

NO IMAGE
AVAILABLE

P26404 : BAYONET SEPARATOR

Curved blade, 1.5mm, 6 3/4" (17cm)

NO IMAGE
AVAILABLE

P26405 : BAYONET SEPARATOR

Curved blade, 1.25mm, 6 3/4" (17cm)

NO IMAGE
AVAILABLE

P26406 : "GIMMICK" ELEVATOR

Angled 15 deg. 1 x 12mm tip

NO IMAGE
AVAILABLE

P26407 : CURETTE

Double-ended, lightly angled, A-1.0mm / B-1.0mm cups

NO IMAGE
AVAILABLE

P26408 : CURETTE

Double-ended, strongly angled, A-1.0mm / B-1.0mm cups

NO IMAGE
AVAILABLE

P26409 : CURETTE

Double-ended, lightly angled, A-2.0mm / B-1.7mm cups

NO IMAGE
AVAILABLE

P26410 : CURETTE

Double-ended, strongly angled, A-2.0mm / B-2.0mm cups

NO IMAGE
AVAILABLE

P26411 : STAPES CURETTE

Double-ended, sharp/blunt, A-2.0mm / B-1.7mm, 6" (15cm)



P26412 : SHORT, STRONG NEEDLE

Curved, 6 1/2" (16.5cm)



P26413 : LIGHT NEEDLE

Curved, 6 1/2" (16.5cm)




P26416 : BARBARA PICK

Straight point, 6 1/2" (16.5cm)




P26417 : OVAL WINDOW PICK

Angled 30 deg., 0.33mm point, 6 1/2" (16.5cm)




P26418 : OVAL WINDOW PICK

Angled 90 deg., fine, 0.33mm point



P26419 : OVAL WINDOW PICK

Angled 90 deg., medium 0.33mm point



P26422 : OVAL WINDOW PICK

Long, angled 90 deg., 0.66mm point, 6 1/2" (16.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

P26423 : HOUGH EXCAVATOR

Angled 90 deg., 0.33 long tip, 6-1/2" (16.5 cm)

NO IMAGE
AVAILABLE

P26424 : PICK, MEDIUM

Angled 90 deg., 1mm point, 6 3/8" (16cm)

NO IMAGE
AVAILABLE

P26425 : PICK, MEDIUM

Angled 90 deg., 1.5mm point, 6 3/8" (16cm)

NO IMAGE
AVAILABLE

P26426 : PICK, MEDIUM

Angled 90 deg., 3mm point, 6 5/16" (16cm)

NO IMAGE
AVAILABLE

P26429 : CRABTREE DISSECTING PICK, ORIGINAL

2.5mm tip, 6 3/8" (16cm)