

NO IMAGE
AVAILABLE

P28356 : Goodyear Uvula Retractor

15mm wide, 7 3/8" (18.5cm)

NO IMAGE
AVAILABLE

P28403 : SNARE WIRE, SIZE 7

one dozen

NO IMAGE
AVAILABLE

P28411 : TONSIL NEEDLE

Straight, 23 Ga., needle extension, 19mm

NO IMAGE
AVAILABLE

P28630 : FRAENKEL LARYNGEAL FORCEPS

Lateral, 8"

NO IMAGE
AVAILABLE

P28774 : Freeman Rake (Yoncoskie) Retractor

7-prong, 3" (7.5cm)

NO IMAGE
AVAILABLE

P28832 : ROSEN CARTILAGE KNIFE

Double-ended, blades set at 45 and 90 degrees, 7 5/8" (19.5cm)

NO IMAGE
AVAILABLE

P28968 : FOX DERMAL CURETTE

Round heads, 1mm dia., 5 1/2" (14cm)

NO IMAGE
AVAILABLE

P28971 : FOX DERMAL CURETTE, #2

Round heads, 4mm dia., 5 1/2" (14cm)

NO IMAGE
AVAILABLE

P28972 : FOX DERMAL CURETTE

Round heads, 5mm dia., 5 1/2" (14cm)

NO IMAGE
AVAILABLE

P28986 : SAALFIELD COMEDONE EXTRACTOR

Double-ended, 5 3/4" (14.5cm)

NO IMAGE
AVAILABLE

P33118 : PIFFARD DERMAL CURETTE, #3

6mm dia., 5 1/2" (14cm)

NO IMAGE
AVAILABLE

P34354 : COMEDONE EXTRACTOR

NO IMAGE
AVAILABLE

P356PRA : Pilling® Blue Fiber Optic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. P356PRA: Pilling Blue Fiber Optic Cable, Pilling illuminator end, 3.5 mm bundle, 6 ft. length, Pilling rotating instrument end

NO IMAGE
AVAILABLE

P358PRA : Pilling® Blue Fiber Optic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes.P358PRA: Pilling Blue Fiber Optic Cable, Pilling illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end

NO IMAGE
AVAILABLE

P88097 : ENT Double Decker Tray with Mat and Inserts

11" x 23" x 5"

NO IMAGE
AVAILABLE

S358PRA : Pilling® Blue Fiber Optic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes.S358PRA: Pilling Blue Fiber Optic Cable, Storz illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end

NO IMAGE
AVAILABLE

W3510LWP : Headlight Cable, 10', 3.5 mm, Bifurcated Lensed (Wolf)

NO IMAGE
AVAILABLE

W3510XLGP : Headlight Cable, 10', 3.5 mm, XLG Lensed (Wolf)

NO IMAGE
AVAILABLE

W356PXA : Pilling® Bifurcated Endfitting Blue Fiberoptic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes.W356PXA: Pilling Blue Fiber Optic Cable, Bifurcated cable, Wolf illuminator end, 3.5 mm bundle, 6 ft. length, Pilling rotating instrument end

NO IMAGE
AVAILABLE

W358LWP : Headlight Cable, 8', 3.5 mm, Bifurcated Lensed (Wolf)

AVAILABLE

NO IMAGE
AVAILABLE

W358PRA : Pilling® Blue Fiber Optic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. W358PRA: Pilling Blue Fiber Optic Cable, Wolf illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end

NO IMAGE
AVAILABLE

W358PXA : Pilling® Bifurcated Endfitting Blue Fiberoptic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. W358PXA: Pilling Blue Fiber Optic Cable, Bifurcated cable, Wolf illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end

NO IMAGE
AVAILABLE

W358XLGP : Headlight Cable, 8', 3.5 mm, XLG Lensed (Wolf)