

NO IMAGE
AVAILABLE

431110 : VIENNA NASAL SPECULUM, SMALL

5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

507485S : MICRO ATOMIZER (SM #'S 151 & 1-RD)

NO IMAGE
AVAILABLE

507611 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

52 Fr., Rounded

NO IMAGE
AVAILABLE

P27814 : SUCTION TUBE, SHORT

3mm, curved

NO IMAGE
AVAILABLE

507615 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

60 Fr., Rounded

NO IMAGE
AVAILABLE

**518446 : JACKSON IMPROVED TRACHEOSTOMY TUBE #8,
W/15MM ADAPTER, REGULAR LENGTH**

9.2mm I.D., 12mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

507608 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

46 Fr., Rounded

NO IMAGE
AVAILABLE

518055 : Jackson Improved Tracheostomy Tube #10

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518055 (size 10, 14.0 mm O.D., 10.6 mm...

NO IMAGE
AVAILABLE

P26769 : SHEA EAR SPECULUM HOLDER

11 3/4" (30cm)

NO IMAGE
AVAILABLE

20121 : ALLIGATOR FORCEPS

Serrated jaws, 4.5 mm, 3-1/8" (8 cm)

NO IMAGE
AVAILABLE

502574 : TYGON TUBING 4MM W/FITTING

NO IMAGE
AVAILABLE

507612 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

54 Fr., Rounded

NO IMAGE
AVAILABLE

518035A : JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #6 FOR SHIKANI-FRENCH SPEAKING VALVE

7.2mm I.D., 10.0mm O.D., 69mm

NO IMAGE
AVAILABLE

507876 : TUCKER LUMEN FINDER

16 Fr., 30 3/4" (78cm)

NO IMAGE
AVAILABLE

518140 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #7

8.1mm I.D., 11mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

025169 : FARRIER OVAL SPECULUM

7.5 x 8.5 mm

NO IMAGE
AVAILABLE

507614 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

58 Fr., Rounded

NO IMAGE
AVAILABLE

518020 : Jackson Improved Tracheostomy Tube #3

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518020 (size 3, 7.0 mm O.D., 4.3 mm...

NO IMAGE
AVAILABLE

507335 : BONNIER CATHETER ADAPTER

1 3/8" (3.5cm)

NO IMAGE
AVAILABLE

507600 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

30 Fr., Rounded

AVAILABLE

NO IMAGE
AVAILABLE

**518441 : JACKSON IMPROVED TRACHEOSTOMY TUBE #7,
W/15MM ADAPTER, REGULAR LENGTH**

8.1mm I.D., 11mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

029480 : ENDAURAL EAR SPECULUM

26mm blade, Strongly curved shaft, 5 5/8" (14.5cm)

NO IMAGE
AVAILABLE

507607 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

44 Fr., Rounded

NO IMAGE
AVAILABLE

509216 : OPERATING PLATFORM

Angled, Ebonized, 0.3" (8mm), 11 5/8" (29.5cm)

NO IMAGE
AVAILABLE

519010 : JACKSON LARYNGECTOMY TUBE #8

10.0mm I.D., 13.6mm O.D., 55mm, stainless steel