

NO IMAGE
AVAILABLE

507877 : TUCKER LUMEN FINDER

18 Fr., 30 3/4" (78cm)

NO IMAGE
AVAILABLE

509207 : VOCAL CORD ROLLERS

Ebonized

NO IMAGE
AVAILABLE

**518904 : JACKSON TRACHEOSTOMY TUBE #7 W/ROTATING LOCK,
EXTRA LONG**

8.3mm I.D., 11mm O.D., 85mm, stainless steel

NO IMAGE
AVAILABLE

507875 : TUCKER LUMEN FINDER

14 Fr., 30 3/4" (78cm)

NO IMAGE
AVAILABLE

511729 : TUCKER VALVE INNER TUBE TO FIT #4 TUCKER TUBE

5.3mm I.D., 58mm, sterling silver

NO IMAGE
AVAILABLE

519030 : JACKSON LARYNGECTOMY TUBE #12

12.7mm I.D., 16.4mm, 55mm, stainless steel

NO IMAGE
AVAILABLE

431112 : VIENNA NASAL SPECULUM, LARGE

5 3/8" (13.5cm)

NO IMAGE
AVAILABLE

506479 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...

NO IMAGE
AVAILABLE

507599 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

28 Fr., Rounded

NO IMAGE
AVAILABLE

514430 : Bernay Tracheal Retractor

Sharp, single prong, 9.5mm depth blade, 6 5/8" (17.0cm)

NO IMAGE
AVAILABLE

P26511 : PICK

3mm curve from tip, 6 5/16" (16cm)

NO IMAGE
AVAILABLE

506802 : Feder-Ossoff Micro-Laryngeal Scissors

Straight, Vertical, Blunt, 9" (23.0cm)

NO IMAGE
AVAILABLE

518030A : JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #5 FOR SHIKANI-FRENCH SPEAKING VALVE

6.2mm I.D., 9.0mm O.D., 68mm

NO IMAGE
AVAILABLE

518240 : JACKSON IMPROVED TRACHEOSTOMY TUBE #8, SHORT LENGTH

9.2mm I.D., 12mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

P27718 : COTTLE SPECULUM

35mm blades, tapered 10mm to 8mm

NO IMAGE
AVAILABLE

507597 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

24 Fr., Rounded

NO IMAGE
AVAILABLE

518139 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #6

7.2mm I.D., 10mm O.D., 69mm, fenestrated, inner and outer

NO IMAGE
AVAILABLE

P27721 : COTTLE SPECULUM

70mm blades, tapered 10mm to 8mm

NO IMAGE
AVAILABLE

030100 : TIECK NASAL SPECULUM, INFANT

Blade 18mm, spread 28mm, 5 1/2" (14.0cm)

NO IMAGE
AVAILABLE

411166 : FARRELL APPLICATOR TRIANGULAR

6 3/8" (16.0cm)

AVAILABLE

NO IMAGE
AVAILABLE

507595 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

20 Fr., Rounded

NO IMAGE
AVAILABLE

**518531 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #5, W/15MM
ADAPTER, REGULAR LENGTH**

6.2mm I.D., 9mm O.D., 68mm, stainless steel

NO IMAGE
AVAILABLE

506952 : MICRO LARYNGEAL FORCEPS

Angled left, fenestrated, serrated, 3 mm tip, 9-1/4" (24 cm)

NO IMAGE
AVAILABLE

507596 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

22 Fr., Rounded

NO IMAGE
AVAILABLE

519015 : JACKSON LARYNGECTOMY TUBE #9

10.7mm I.D., 14mm O.D., 55mm, stainless steel