

NO IMAGE
AVAILABLE

506950 : MICRO LARYNGEAL FORCEPS

Angled right, fenestrated, serrated, 3 mm tip, 9-1/4" (24 cm)

NO IMAGE
AVAILABLE

507598 : HURST TUNGSTEN FILLED ESOPHAGEAL BOUGIE

26 Fr., Rounded

NO IMAGE
AVAILABLE

518906 : JACKSON TRACHEOSTOMY TUBE #8 W/ROTATING LOCK, EXTRA LONG

9.2mm I.D., 12mm O.D., 85mm, stainless steel

NO IMAGE
AVAILABLE

514380 : LABORDE DILATOR

5 1/2" (14cm)

NO IMAGE
AVAILABLE

518120 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #3

4.3mm I.D., 7mm O.D., 58mm, stainless steel

NO IMAGE
AVAILABLE

522590 : BRONCHOSCOPE ADAPTER 11MM TO 15MM

3/4" (2.0cm)

NO IMAGE
AVAILABLE

025103 : BOUCHERON EAR SPECULUM, ROUND, SIZE 3

Diameter 4.6mm

NO IMAGE
AVAILABLE

506856 : FEDER-OSSOFF MICRO LARYNGEAL GRASPING FORCEPS

Double action, straight, with suction. Shaft tapers from 2.0 mm to 1 mm. Coarsely serrated jaws. 1.0 mm suction channel w/ luer-lok proximal fitting.

NO IMAGE
AVAILABLE

516390 : MARTIN LARYNGECTOMY TUBE #10

10.6mm I.D., 14mm O.D., 42mm, sterling silver

NO IMAGE
AVAILABLE

518145 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #8

9.2mm I.D., 12mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

025101 : BOUCHERON EAR SPECULUM, ROUND, SIZE 1

Diameter 3mm

NO IMAGE
AVAILABLE

054015 : TYRELL HOOK

Sharp, 1.3 mm hook, 4-3/4" (12 cm)

NO IMAGE
AVAILABLE

506820 : Feder-Ossoff Micro-Laryngeal Scissors

Angled down, Vertical, Sharp, 9" (23.0cm)

NO IMAGE
AVAILABLE

518049 : Jackson Improved Tracheostomy Tube - Fenestrated #8

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...

NO IMAGE
AVAILABLE

P26601 : SHEEHY FASCIA PRESS

23mm x 31mm plate, 7 1/2" (19cm)

NO IMAGE
AVAILABLE

507816 : JACKSON ESOPHAGEAL DILATOR

16 Fr.

NO IMAGE
AVAILABLE

518155 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #10

10.6mm I.D., 14mm O.D., 69mm

NO IMAGE
AVAILABLE

P27720 : COTTLE SPECULUM

55mm blades, tapered 10mm to 8mm

NO IMAGE
AVAILABLE

025104 : BOUCHERON EAR SPECULUM, ROUND, SIZE 4

Diameter 5.1mm

NO IMAGE
AVAILABLE

505110 : JACKSON LARYNGEAL CUP FORCEPS

Straight, 2mm cup, 20 cm

AVAILABLE

NO IMAGE
AVAILABLE

506668 : Jako Sharp Hook, Acute

Dull finish to reduce glare. May be used with 506888 Feder-Ossoff Universal Ring Handled Knife Handle. Working length 26 cm.506668: Jako Sharp Hook, acute

NO IMAGE
AVAILABLE

511139P : TUCKER VALVE INNER TUBE ONLY FOR #6 JACKSON IMPROVED TUBE

7.1mm I.D., 72mm, sterling silver

NO IMAGE
AVAILABLE

514412 : Lukens Thymus Retractor

5/8" (1.5cm) depth blade, 6 7/8" (17.5cm)

NO IMAGE
AVAILABLE

511734 : TUCKER VALVE INNER TUBE TO FIT #5 TUCKER TUBE

6.1mm I.D., 63mm, sterling silver

NO IMAGE
AVAILABLE

511135P : JACKSON IMPROVED TRACHEOSTOMY TUBE WITH ROTATING LOCK #6, INNER AND OUTER TUBES AND OBTURATOR

7.1mm I.D., 10.0mm O.D., 69mm, sterling silver