

NO IMAGE
AVAILABLE

514377 : JACKSON ADULT DILATOR

5 1/2" (14cm)

NO IMAGE
AVAILABLE

518541 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #7, W/15MM ADAPTER, REGULAR LENGTH

8.1mm I.D., 11mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

506667 : Jako Sharp Hook, Slight

Dull finish to reduce glare. May be used with 506888 Feder-Ossoff Universal Ring Handled Knife Handle. Working length 26 cm. 506667: Jako Sharp Hook, slight

NO IMAGE
AVAILABLE

515120 : AIR-LON® TRACHEA TUBE #7

11.0mm O.D.

NO IMAGE
AVAILABLE

025102 : BOUCHERON EAR SPECULUM, ROUND, SIZE 2

Diameter 3.6mm

NO IMAGE
AVAILABLE

026275 : DINGMAN ZYGOMA HOOK

Sharp, 13 mm hook, 7" (17.8 cm)

NO IMAGE
AVAILABLE

506806 : Feder-Ossoff Micro-Laryngeal Scissors

Curved right, Horizontal, Blunt, 9" (23.0cm)

NO IMAGE
AVAILABLE

519000 : JACKSON LARYNGECTOMY TUBE #6

7.2mm I.D., 10mm O.D., 55mm, stainless steel

NO IMAGE
AVAILABLE

524041 : REPLACEMENT TELESCOPE GUIDE RUBBER CAP FOR 524038

For 5.5 mm telescope, 524036, 5 per pack, 5/8" (1.5cm)

NO IMAGE
AVAILABLE

518635 : JACKSON IMPROVED TRACHEOSTOMY TUBE #7, W/15MM ADAPTER, SHORT LENGTH

8.1mm I.D., 11mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

1037833 : DECANNULATION STOPPER

White stopper with 40" cord, Size 3, Full Closure

NO IMAGE
AVAILABLE

351970 : BOVIE ELECTROSURGICAL HANDLE

5 1/8" (13.0cm)

NO IMAGE
AVAILABLE

519025 : JACKSON LARYNGECTOMY TUBE #11

12.3mm I.D., 15.4mm, 55mm, stainless steel

NO IMAGE
AVAILABLE

524040 : REPLACEMENT TELESCOPE GUIDE RUBBER CAP FOR 524037

For 4.0 mm telescope, 524036, 5 per pack, 5/8" (1.5cm)

NO IMAGE
AVAILABLE

507814 : JACKSON ESOPHAGEAL DILATOR

14 Fr.

NO IMAGE
AVAILABLE

518029 : Jackson Improved Tracheostomy Tube - Fenestrated #4

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...

NO IMAGE
AVAILABLE

054016 : TYRELL HOOK

Blunt, 1.3 mm hook, 4-3/4" (12 cm)

NO IMAGE
AVAILABLE

20100 : HOUSE ALLIGATOR FORCEPS

Serrated jaws, 8 mm, 3-1/8" (8 cm)

NO IMAGE
AVAILABLE

507812 : JACKSON ESOPHAGEAL DILATOR

12 Fr.

NO IMAGE
AVAILABLE

511725 : TUCKER TUBE, INNER AND OUTER TUBES, OBTURATOR, SIZE 4

5.3mm I.D., 8mm O.D., 55mm, sterling silver

AVAILABLE

NO IMAGE
AVAILABLE

054720 : Aufricht Nasal Retractor

Small blade 8mm wide, angled shaft

NO IMAGE
AVAILABLE

518325 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #5, SHORT LENGTH

6.2mm I.D., 9mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

507830 : JACKSON ESOPHAGEAL DILATOR

30 Fr.

NO IMAGE
AVAILABLE

518620 : JACKSON IMPROVED TRACHEOSTOMY TUBE #4, W/15MM ADAPTER, SHORT LENGTH

5.5mm I.D., 8mm O.D., 52mm, stainless steel

NO IMAGE
AVAILABLE

507808 : JACKSON ESOPHAGEAL DILATOR

8 Fr.