


General

**443510 : JACOBSON HEMOSTATIC FORCEPS**

Curved, Serrated, 1.2mm tip, 7 1/8" (18 cm)

**510101 : MAYO HEGAR NEEDLE HOLDER**

15.5 x 3mm jaw, 6" (15.5cm), TC

**060100 : Skin Retractor**

Slightly angled, 3 3/4" (9.5cm)

**165382 : BECKMAN-WEITLANER RETRACTOR**


3 x 4 Blunt prongs, 6 1/2" (16.5cm)

**181223 : ADSON TISSUE FORCEPS**

1 x 2 Teeth, 4 3/4" (12 cm)

**510110 : MAYO HEGAR NEEDLE HOLDER**

16.2 x 2.8mm jaw, 7" (18cm), TC



463913 : Ultimate Cut Metzenbaum Scissors

Curved, w/one Serrated blade, 7" (18.0cm)



89-2401 : PITIE-SALPETRIERE VEIN RTR HK 24CM, VEIN HOOK

For use with surgical blades no. 10, 11, 12, and 15.



165358 : WEITLANER RETRACTOR

2 x 3 Blunt prongs, 4 1/2"(11.5cm)



510151 : MAYO HEGAR BROAD BULLDOG NEEDLE HOLDER

11.8 x 4mm jaw, 7" (18cm), TC



806024 : FLASH TRAY

6 1/2" (16.5cm) x 10" (25.5cm)



144300 : Iris Scissors

Straight, two sharp points, 4 1/2" (11.5cm)



054125 : JOHNSON SKIN HOOK

Sharp, 5 mm hook, single prong, 6-1/2" (16.5 cm)

164404 : Volkmann Rake Retractor

Sharp, 4-pronged, 3/4" (2.0cm) width blade, 8 1/2" (21.5cm)



485450 : O'SULLIVAN-O'CONNOR ABDOMINAL RETRACTOR, LARGE

Spread 12 1/2" (31.5cm), Side Blades 1 1/8" (2.7cm) wide x 2 7/8" (73mm) deep, Center blades: two small 1 3/4" (44mm) wide x 3 1/8" (79mm) deep & one large 2 7/8" (72mm) wide x 4" (102mm) deep




165364 : WEITLANER RETRACTOR

3 x 4 Blunt prongs, 7 5/8" (19.5cm)



510420 : CRILE WOOD NEEDLE HOLDER

15.4 x 2.3mm jaw, 10 1/4" (26cm), TC




730295 : Hill Ferguson Retractor

1" (25mm) width x 2 3/8" (6.0cm) depth blade, 8 1/2" (21.5cm)



141312 : Mayo Scissors

Straight, 6 3/4" (17.0cm)



202623 : BARR CRYPT HOOK


Double-ended, 10 1/8" (25.65cm)





481230 : Cushing Vain / Nerve Retractor

5/8" (16mm) width x 1/2" (1.5cm) depth blade, 8 5/8" (22.0cm)



060101 : Skin Retractor

Strongly angled, 3 1/2" (9.0cm)



KM52486 : Beckman-Adson Laminectomy Retractor

4 x 4 prongs, hinged blades 25 mm high, 12-1/4" (31 cm). KM52486: Sharp, 44 mm deep



508004 : Video Laryngoscope Tray with Mat

Pilling Trays are designed for microsurgical and other delicate instruments. These lightweight, durable trays are engineered for maximum protection and flexibility. The transparent lid will help you easily identify...



165360 : WEITLANER RETRACTOR

3 x 4 Blunt prongs, 5 1/2" (14cm)