

General



NO IMAGE
AVAILABLE

142220 : Mayo Scissors

Curved, 5 1/2" (14.0cm), TC



NO IMAGE
AVAILABLE

164166 : Ochsner Ribbon Retractor

1 1/4" (3.2cm) width blade, 13" (33.0cm)



NO IMAGE
AVAILABLE

463941 : Ultimate Cut Potts Scissors

Angled 45, w/one Serrated blade, 7 1/4" (18.5cm)



NO IMAGE
AVAILABLE

P13082 : SELMAN TISSUE FORCEPS

Round cup w/ teeth, 12" (30.5 cm)



NO IMAGE
AVAILABLE

462641 : Iris Scissors

Curved, 5" (12.5cm)



NO IMAGE
AVAILABLE

476330 : CARMALT FORCEPS

Curved, 3mm tip, 6 1/2" (17 cm)

NO IMAGE
AVAILABLE

P13167 : BAYONET NEEDLE HOLDER

14mm jaws, 7 1/4" (18.5cm)

NO IMAGE
AVAILABLE

480172 : Parker Retractor, Large

Double ended, 1" (2.5cm) / 1" (2.5cm) width blade, 7 1/8" (18.0cm), 2/set

NO IMAGE
AVAILABLE

P16532 : DIETRICH NEEDLE HOLDER

Serrated jaws, 7" (18cm), TC

NO IMAGE
AVAILABLE

P16534 : DIETRICH NEEDLE HOLDER

Serrated jaws, 9" (23cm), TC

NO IMAGE
AVAILABLE

454100 : FOERSTER SPONGE FORCEPS

Straight, Serrated, 9 1/2" (24 cm)

NO IMAGE
AVAILABLE

P16544 : JACOBSON NEEDLE HOLDER, HEAVY

Round handle, w/lock, 7 1/2" (19cm), TC

NO IMAGE
AVAILABLE

P16545 : JACOBSON NEEDLE HOLDER, HEAVY

w/lock, 9" (23cm), TC

NO IMAGE
AVAILABLE

164456 : Murphy Retractor

Blunt, 6-pronged, 1 1/8" (3.0cm) width blade, 7 3/4" (19.5cm)

NO IMAGE
AVAILABLE

464502 : Operating Scissors

Straight, Sharp/Sharp blades, 5 1/2" (14.0cm), TC

NO IMAGE
AVAILABLE

P16551 : JACOBSEN MICRO NEEDLE HOLDER

w/lock, 7 1/4" (18.5cm), TC

NO IMAGE
AVAILABLE

460841 : Metzenbaum Scissors

Weck Pattern, bright finish. 460841: Curved, 5-3/8" (13.5 cm)

NO IMAGE
AVAILABLE

P16552 : JACOBSEN MICRO NEEDLE HOLDER

w/lock, 8 1/4" (21cm), TC

NO IMAGE
AVAILABLE

P16553 : JACOBSEN MICRO NEEDLE HOLDER

Round handle, w/lock, 9" (23cm), TC

NO IMAGE

P16555 : JACOBSEN MICRO NEEDLE HOLDER

straight, w/o lock, 9" (23cm), TC

AVAILABLE

NO IMAGE
AVAILABLE

464721 : Metzenbaum Scissors

Straight, 9" (23.0cm), Delicate, TC

NO IMAGE
AVAILABLE

P16734 : MAYO HEGAR NEEDLE HOLDER, LEFT-HANDED

10" (25.5cm), TC

NO IMAGE
AVAILABLE

476220 : OCHSNER FORCEPS

Straight, Full serrations, 1 x 2 teeth, 72 x 5.5mm tip, 8 7/8" (23 cm)

NO IMAGE
AVAILABLE

P16740 : DEBAKEY ULTRA LIGHT NEEDLE HOLDER

7" (18cm), TC

NO IMAGE
AVAILABLE

462691 : Strabismus Scissors

Curved, 4 1/2" (11.5cm)