

## General



NO IMAGE  
AVAILABLE

**200234 : ROSEBUDS (COTTON), 3/8"**

X-Ray Detectable, Box of 40 Packages, Sterile, Cotton, 5/PKG



NO IMAGE  
AVAILABLE

**200235 : ROSEBUDS (COTTON), 1/4"**

X-Ray Detectable, Box of 20 Packages, Sterile, Strung, Gauze, 10/PKG



NO IMAGE  
AVAILABLE

**200236 : ROSEBUDS (COTTON), 1/2"**

X-Ray Detectable, Box of 20 Packages, Sterile, Strung, Gauze, 10/PKG



NO IMAGE  
AVAILABLE

**200240 : BRONCHOSCOPIC SPONGES (COTTON), 1/4" x 3/4"**

X-Ray Detectable, Box of 20 Packages, Sterile, Cotton, 10 Sponges on Safety Pin per Package



NO IMAGE  
AVAILABLE

**200250 : SPONGES (COTTON), 2" x 2"**

X-Ray Detectable, Box of 20 Packages, Sterile, Strung, Gauze, 10/PKG



NO IMAGE  
AVAILABLE

**200322 : TONSIL SPONGES (COTTON), LARGE**

X-Ray Detectable, Qty 5, Non-sterile, Double-strung, Gauze, 1000/PKG

NO IMAGE  
AVAILABLE

**200326 : TONSIL SPONGES (COTTON), MEDIUM**

X-Ray Detectable, Box of 100 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200327 : TONSIL SPONGES (COTTON), LARGE**

X-Ray Detectable, Box of 100 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200328 : TONSIL SPONGES (COTTON), X-LARGE**

X-Ray Detectable, Box of 100 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200333 : TONSIL SPONGES, MEDIUM**

X-Ray Detectable, Box of 100 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200334 : TONSIL SPONGES, LARGE**

X-Ray Detectable, Box of 100 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200335 : TONSIL SPONGES, X-LARGE**

X-Ray Detectable, Box of 100 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200340 : TONSIL SPONGES, MEDIUM**

X-Ray Detectable, Box of 20 Packages, Sterile, Strung, Gauze, 5/PKG

NO IMAGE  
AVAILABLE

**200400 : CYLINDRICAL SPONGES, 3/8" x 1 1/2"**

Non X-Ray Detectable, Box of 20 Packages, Sterile, End-strung, Cotton, 5/PKG

NO IMAGE  
AVAILABLE

**200401 : CYLINDRICAL SPONGES, 3/8" x 1 1/2"**

Non X-Ray Detectable, Box of 20 Packages, Sterile, Center-strung, Cotton, 5/PKG

NO IMAGE  
AVAILABLE

**200402 : CYLINDRICAL SPONGES, 1/2" x 1 1/2"**

Non X-Ray Detectable, Box of 20 Packages, Sterile, End-strung, Cotton, 5/PKG

NO IMAGE  
AVAILABLE

**200403 : CYLINDRICAL SPONGES, 1/2" x 1 1/2"**

Non X-Ray Detectable, Box of 20 Packages, Sterile, Center-strung, Cotton, 5/PKG

NO IMAGE  
AVAILABLE

**200550 : WECK® Wipe**

X-ray detectable, sterile. 20 per box, 1 per pouch

NO IMAGE  
AVAILABLE

**202030 : BUIE CLAMP**

Curved, Serrated, 5/8" (1.6cm) tip length, 8 7/8" (22.5cm)

NO IMAGE  
AVAILABLE

**202095 : HIRSCHMAN ANOSCOPE**

Small, 5/8" (16mm) tip diameter, 2 1/2" (6.5cm)

AVAILABLE

NO IMAGE  
AVAILABLE

**202096 : HIRSCHMAN ANOSCOPE**

Medium, 3/4" (19mm) tip diameter, 2 3/4" (7cm)

NO IMAGE  
AVAILABLE

**202097 : HIRSCHMAN ANOSCOPE**

Large, 7/8" (22mm) tip diameter, 2 3/4" (7cm)

NO IMAGE  
AVAILABLE

**202600 : LARRY PROBE AND GROOVE DIRECTOR**

Stainless Steel, 5 7/8" (15cm)

NO IMAGE  
AVAILABLE

**202630 : STEWART RECTAL HOOK**

7/8" (22 mm) hook, 8-7/8" (22.7 cm)

NO IMAGE  
AVAILABLE

**202792 : BUIE SUCTION TUBE w/ THUMB VALVE**

15 Fr. Diameter, 17" (43.2 cm)