

General



534201 : PATTERSON TROCAR SIZE 1

7 3/4" (19.5cm)




534202 : PATTERSON TROCAR SIZE 2

8" (20.5cm)




534203 : PATTERSON TROCAR SIZE 3

8" (20.5cm)




536100 : ADSON SUCTION TUBE

Curved, 14 French tip, 7 1/4" (18.5cm)



536103 : ADSON SUCTION TUBE

Straight, 12 French tip, 8 1/4" (21.0cm)



536104 : ADSON SUCTION TUBE

Curved, 8 1/4" (21.0cm)

NO IMAGE
AVAILABLE

536105 : ADSON SUCTION TUBE

Curved, 6 1/2" (16.5cm)

NO IMAGE
AVAILABLE

536107 : POOLE PEDIATRIC SUCTION

16 French, 6mm tip, Working length 7 1/2" (19cm)

NO IMAGE
AVAILABLE

536220 : ANDREWS PYNCHON SUCTION TUBE

Curved, 6mm tip, 3mm dia. removable shaft, 9 5/8" (24.5cm)

NO IMAGE
AVAILABLE

550190 : MASTIN MUSCLE FORCEPS

Curved on side, 3.1 mm tip, 6-5/8" (17 cm)

NO IMAGE
AVAILABLE

65216 : KLEINERT-KUTZ DUAL ACTION BONE CUTTING FORCEPS

Straight

NO IMAGE
AVAILABLE

658505 : Raney Nylon Clips

[href="http://www.teleflexsurgicalcatalog.com/surgical/pilling/product/658505"](http://www.teleflexsurgicalcatalog.com/surgical/pilling/product/658505)658505:

Sterile, single use, 10/pouch, 20 pouches/box

NO IMAGE
AVAILABLE

746550 : HERRICK KIDNEY CLAMP

Angled, Longitudinal Serrations, 1 5/8" (4.1cm) tip length, 9" (23cm)

NO IMAGE
AVAILABLE

746797 : GUYON CATHETER GUIDE

Curved, 9 Fr., 13 1/4" (33.5cm)

NO IMAGE
AVAILABLE

746798 : VAN BUREN CATHETER GUIDE

Curved, 9 Fr., 13 3/4" (35.0cm)

NO IMAGE
AVAILABLE

746800 : ZIPSER PENIS CLAMP

3 3/8" (8.5cm)

NO IMAGE
AVAILABLE

89-2712 : SAPHLITE BLADE, LARGE

NO IMAGE
AVAILABLE

89-2713 : SAPHLITE BLADE, MID

NO IMAGE
AVAILABLE

89-2714 : SAPHLITE BLADE, MINI

NO IMAGE

89-2715 : SAPHLITE BLADE, SKINI-MINI

AVAILABLE

NO IMAGE
AVAILABLE

89-2733 : SAPHLITE BLADE, SKINI-MID

NO IMAGE
AVAILABLE

89-2734K : SAPHLITE BLADE, COBRA

NO IMAGE
AVAILABLE

89-2750K : SAPHLITE BLADE, CURVED MID

NO IMAGE
AVAILABLE

89-9463K : BALL JOINT FOR GENZARM SYSTEM

NO IMAGE
AVAILABLE

89-9489 : GENZARM BALL TIGHTENER