

General



KM38094 : DRESSING FORCEPS

8"




KM38100 : DRESSING FORCEPS

10"




KM38106 : DRESSING FORCEPS

12"




KM38180 : GERALD DRESSING FORCEPS

Straight, serrated, 7"



KM38184 : CUSHING BAYONET DRESSING FORCEPS

Serrated, 7-1/4" (18.5 cm)



KM38186 : CUSHING BAYONET TISSUE FORCEPS

1 x 2 teeth, 7-1/4"

NO IMAGE
AVAILABLE

KM38264 : TISSUE FORCEPS

1 x 2, 4-1/2"

NO IMAGE
AVAILABLE

KM38268 : TISSUE FORCEPS

1 x 2, 5"

NO IMAGE
AVAILABLE

KM38272 : TISSUE FORCEPS

1 x 2, 5-1/2"

NO IMAGE
AVAILABLE

KM38276 : TISSUE FORCEPS

1 x 2, 6"

NO IMAGE
AVAILABLE

KM38280 : TISSUE FORCEPS

1 x 2, 8"

NO IMAGE
AVAILABLE

KM38286 : TISSUE FORCEPS

1 x 2, 10"

NO IMAGE
AVAILABLE

KM38290 : TISSUE FORCEPS

1 x 2, 12"

NO IMAGE
AVAILABLE

KM38296 : THUMB TISSUE FORCEPS

2 x 3 teeth, 3.3 mm tip, 5-1/2" (14 cm)

NO IMAGE
AVAILABLE

KM38322 : FERRIS SMITH TISSUE FORCEPS

1 x 2 teeth, 7"

NO IMAGE
AVAILABLE

KM38324 : FERRIS-SMITH TISSUE FORCEPS

Serrated, 2 x 3 teeth, 3.9 mm tip, 7" (18 cm)

NO IMAGE
AVAILABLE

KM39052 : ADSON DRESSING FORCEPS

Serrated, 1.2 mm tip, 6" (15 cm)

NO IMAGE
AVAILABLE

KM39063 : ADSON TISSUE FORCEPS

Serrated, 1 x 2 teeth, 4-3/4" (12 cm)

NO IMAGE
AVAILABLE

KM39066 : ADSON FORCEPS

1 x 2 standard hudson, 4-3/4"

NO IMAGE
AVAILABLE

KM39068 : ADSON FORCEPS

2 x 3, 4-3/4"

AVAILABLE

NO IMAGE
AVAILABLE

KM39070 : ADSON TISSUE / SUTURE FORCEPS

1 x 2, tying platform

NO IMAGE
AVAILABLE

KM39072 : ADSON FORCEPS

Serrated, carbide, 4-3/4"

NO IMAGE
AVAILABLE

KM39074 : ADSON BAYONET FORCEPS

Straight, serrated, 8-1/4"

NO IMAGE
AVAILABLE

KM39076 : SEMKEN DRESSING FORCEPS

Serrated, 5"

NO IMAGE
AVAILABLE

KM39080 : SEMKEN DRESSING FORCEPS

Straight, 1 x 2 teeth, 1.2 mm tip, 5" (13 cm)