

Home // General

General

NO IMAGE
AVAILABLE

KM41175L : CRILE-WOOD NEEDLE HOLDER

Left-hand, Serrated, 7", TC

NO IMAGE
AVAILABLE

KM41176 : CRILE-WOOD NEEDLE HOLDER

Serrated, 8", TC

NO IMAGE
AVAILABLE

KM41178 : CRILE-WOOD NEEDLE HOLDER

Serrated, 9", TC

NO IMAGE
AVAILABLE

KM41257 : RYDER NEEDLE HOLDER

5", TC

NO IMAGE
AVAILABLE

KM41264 : RYDER NEEDLE HOLDER

6", TC

NO IMAGE
AVAILABLE

KM41265 : RYDER NEEDLE HOLDER

7", TC

NO IMAGE
AVAILABLE

KM41266 : RYDER NEEDLE HOLDER

8", TC

NO IMAGE
AVAILABLE

KM41268 : RYDER NEEDLE HOLDER

10 1/4", TC

NO IMAGE
AVAILABLE

KM41302 : MAYO HEGAR NEEDLE HOLDER

Serrated, 6", TC

NO IMAGE
AVAILABLE

KM41302L : MAYO HEGAR NEEDLE HOLDER

Left-hand, Serrated, 6", TC

NO IMAGE
AVAILABLE

KM41308 : MAYO HEGAR NEEDLE HOLDER

Serrated, 7", TC

NO IMAGE
AVAILABLE

KM41314 : MAYO HEGAR NEEDLE HOLDER

8", TC

NO IMAGE
AVAILABLE

KM41316 : MAYO-HEGAR (MASSON) NEEDLE HOLDER

10 1/2", TC

NO IMAGE
AVAILABLE

KM41333 : GILLIES NEEDLE HOLDER

Curved, 6 1/2"

NO IMAGE
AVAILABLE

KM41448 : OLSEN-HEGAR NEEDLE HOLDER

5 1/2", TC

NO IMAGE
AVAILABLE

KM41454 : OLSEN-HEGAR NEEDLE HOLDER

7 1/2", TC

NO IMAGE
AVAILABLE

KM42132 : HEANEY NEEDLE HOLDER

8 1/4"

NO IMAGE
AVAILABLE

KM42134 : HEANEY NEEDLE HOLDER

Curved jaw, 8 1/4", TC

NO IMAGE
AVAILABLE

KM42152 : SAROT NEEDLE HOLDER

7 1/2", TC

NO IMAGE
AVAILABLE

KM42254 : MICHEL CLIP APPLYING FORCEPS

5" (13 cm)

AVAILABLE

NO IMAGE
AVAILABLE

**KM42256 : WACHENFELDT CLIP APPLYING AND REMOVING
FORCEPS**

4.1 mm tip, 5" (13 cm)

NO IMAGE
AVAILABLE

**KM42258 : WACHENFELDT CLIP APPLYING AND REMOVING
FORCEPS**

3.5 mm tip, 5" (13 cm)

NO IMAGE
AVAILABLE

KM42608 : PROBE WITH EYE

Malleable, stainless, 4 1/2"

NO IMAGE
AVAILABLE

KM42610 : PROBE WITH EYE

Malleable, stainless, 5"

NO IMAGE
AVAILABLE

KM42612 : PROBE WITH EYE

Malleable, stainless, 5 1/2"