

## General

NO IMAGE  
AVAILABLE

**KM53380 : IRIS FORCEPS 4"**

1 x 2, straight

NO IMAGE  
AVAILABLE

**KM53450 : CASTRO SUTURE FORCEPS**

Straight, 1 x 2 teeth, 0.12 mm, 4"

NO IMAGE  
AVAILABLE

**KM57164 : GRUENWALD (JANSEN) BAYONET FORCEPS**

8"

NO IMAGE  
AVAILABLE

**KM57178 : WILDE EAR FORCEPS**

Serrated, 5" (13 cm)

NO IMAGE  
AVAILABLE

**KM58500 : MAGILL CATHETER FORCEPS, CHILD**

7"

NO IMAGE  
AVAILABLE

**KM58502 : MAGILL CATHETER FORCEPS, ADULT**

9"

NO IMAGE  
AVAILABLE

**KM59064 : SCHNIDT FORCEPS**

Half-curved, closed ring, 7-1/2"

NO IMAGE  
AVAILABLE

**KM59066 : SCHNIDT FORCEPS**

Full-curved, closed ring, 7-1/2"

NO IMAGE  
AVAILABLE

**KM63300 : BODENHAMMER REC.SPECULUM 6" LONG**

NO IMAGE  
AVAILABLE

**KM63354 : BRINKERHOFF RECTAL SPECULUM, LARGE**

NO IMAGE  
AVAILABLE

**KM63503 : MCGIVNEY FORCEPS**

7-1/2" (19 cm)

NO IMAGE  
AVAILABLE

**KM66058 : BOZEMAN UTERINE PACKING FORCEPS**

Curved, 10"

NO IMAGE  
AVAILABLE

**KM67500 : DE LEE OBSTETRICAL FORCEPS**

Fenestrated blades, long model, 14" (36 cm)

NO IMAGE  
AVAILABLE

**KM67502 : SIMPSON OBSTETRIC FORCEPS**

12"

NO IMAGE  
AVAILABLE

**KM83512 : HALSTEAD TITANIUM MOSQUITO FORCEPS**

Curved, 5" (13 cm)

NO IMAGE  
AVAILABLE

**KM83522 : CRILE TITANIUM HEMOSTAT FORCEPS**

Curved, 2.1 mm tip, 5-1/2" (14 cm)

NO IMAGE  
AVAILABLE

**KM83544 : ALLIS TITANIUM TISSUE FORCEPS 4X5**

4 x 5 teeth, 6" (15 cm)

NO IMAGE  
AVAILABLE

**KM83547 : BACKHAUS TOWEL FORCEPS**

Curved, 4-1/8" (11 cm)

NO IMAGE  
AVAILABLE

**P12481 : KNIFE HANDLE**

Round, knurled, chuck style

NO IMAGE  
AVAILABLE

**P12483 : KNIFE HANDLE, LONG**

Round, knurled, chuck style, 6" (15cm)

AVAILABLE

NO IMAGE  
AVAILABLE

**P12840 : OCHSNER TROCAR, 16 FR**

For 12 Fr catheter

NO IMAGE  
AVAILABLE

**P12842 : OCHSNER TROCAR, 20 FR**

For 16 Fr catheter

NO IMAGE  
AVAILABLE

**P12857 : O'GRADY CORONARY ARTERY SUCTION TIP**

12 1/2" (32cm)

NO IMAGE  
AVAILABLE

**P12868 : DEBAKEY-ADSON SUCTION TUBE**

4" working length, 6 1/16" (15.5cm)

NO IMAGE  
AVAILABLE

**P12869 : DEBAKEY-ADSON SUCTION TUBE**

6" working length, 8 1/16" (20.5cm)