

General

NO IMAGE
AVAILABLE**P37333 : Surgical Case**

Case with optional trays and support brackets, without pin mat. P37333 ext. size: 11" x 23" x 5" (28 x 58.5 x 12.5 cm) int. size: 10.7" x 21.62" x 4.91" (27 x 55 x 12.5 cm)

NO IMAGE
AVAILABLE**P37334 : Surgical Case**

Case with 2 optional trays and support brackets, without pin mat. P37334 ext. size: 11" x 23" x 5" (28 x 58.5 x 12.5 cm) int. size: 10.7" x 21.62" x 4.91" (27 x 55 x 12.5 cm)

NO IMAGE
AVAILABLE**P37335 : Osteotome Case for Lambottes**

For curved and straight Lambottes, ext. size: 9" x 11-1/2" x 3" (23 x 29 x 7.5 cm) int. size: 8.57" x 11.12" x 2.38" (22 x 28 x 6 cm)

NO IMAGE
AVAILABLE**P37336 : Surgical Utility Case**

Ext. size: 5" x 10-1/2" x 2" (12.5 x 26.5 x 5 cm) int. size: 4.5" x 10.29" x 1.63" (11.5 x 26 x 4 cm)

NO IMAGE
AVAILABLE**P37337 : Surgical Utility Case**

With full pin mat, center lid retainer, ext. size: 5" x 10-1/2" x 2" (12.5 x 26.5 x 5 cm) int. size: 4.5" x 10.29" x 1.53" (11.5 x 26 x 4 cm)

NO IMAGE
AVAILABLE**P37338 : Surgical Tray**

P37338: 10-1/2" x 21-1/2" x 2-1/4" (26.5 x 54.5 x 5.5 cm)

NO IMAGE
AVAILABLE

P37339 : Universal Keyhole Bracket

10 position, 7-1/2" (19 cm)

NO IMAGE
AVAILABLE

P37341 : Universal Keyhole Bracket

Universal Keyhole Bracket, 9-1/4" (23.5 cm)

NO IMAGE
AVAILABLE

P37344 : Slotted "U" Bracket

10 mm slot, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37345 : Slotted "U" Bracket

12 mm slot, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37349 : Slotted "U" Bracket

10 mm slot, 1" (2.5 cm)

NO IMAGE
AVAILABLE

P37351 : Pass Thru Bracket

6 mm hole, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37352 : Pass Thru Bracket

7 mm hole, 1-1/2" (4 cm)



P37353 : Pass Thru Bracket

10 mm hole, 1-1/2" (4 cm)



P37355 : Pass Thru Bracket

6 mm hole, 1" (2.5 cm)



P37356 : Pass Thru Bracket

10 mm hole, 1" (2.5 cm)



P37357 : Pass Thru Bracket

12 mm hole, 1" (2.5 cm)



P37358 : Universal Keyhole Bracket

6 position, 5" (12.5 cm)



P37372 : Keyhole Bracket

3 mm hole, 1-1/2" (4 cm)



P37373 : Keyhole Bracket

4 mm hole, 1-1/2" (4 cm)

AVAILABLE

NO IMAGE
AVAILABLE

P37374 : Keyhole Bracket

6 mm hole, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37376 : Keyhole Bracket

10 mm hole, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37378 : Keyhole Bracket

13 mm hole, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37379 : Keyhole Bracket

15 mm hole, 1-1/2" (4 cm)

NO IMAGE
AVAILABLE

P37380 : Keyhole Bracket

19 mm hole, 1-1/2" (4 cm)