

Home // General

## General




### **P37382 : Keyhole Bracket**

25 mm hole, 1-1/2" (4 cm)




### **P37383 : Keyhole Bracket**

4 mm hole, 1" (2.5 cm)




### **P37384 : Keyhole Bracket**

5 mm hole, 1" (2.5 cm)




### **P37387 : Keyhole Bracket**

10 mm hole, 1" (2.5 cm)



### **P37390 : Keyhole Bracket**

3 tier 10/5/3 mm, 1-1/2" (4 cm)



### **P37395 : Pin Mat**

9.63" x 14.13"

NO IMAGE  
AVAILABLE

**P37396 : Pin Mat**

9.63" x 18"

NO IMAGE  
AVAILABLE

**P37397 : Pin Mat**

9.5" x 9.5"

NO IMAGE  
AVAILABLE

**P37398 : Pin Mat**

10.12" x 18"

NO IMAGE  
AVAILABLE

**P37404 : Pin Mat**

13.13" x 18.25"

NO IMAGE  
AVAILABLE

**P37419 : Surgical Case with Optional Tray**

Case without brackets/pin mat. P37419 ext. size: 11" x 23" x 3-1/2" (28 x 58.5 x 9 cm) int. size: 10.70" x 21.62" x 3.41" (27.5 x 55 x 8.5 cm)

NO IMAGE  
AVAILABLE

**P37432 : Keyhole Bracket**

2 mm hole, 1" (2.5 cm)

NO IMAGE  
AVAILABLE

**P37433 : Keyhole Bracket**

3 mm hole, 1" (2.5 cm)




**P37437 : Pin Mat**

1" x 4.5"




**P37454 : Pin Mat**

4.25" x 9.75"




**P37485 : Keyhole Bracket**

5 mm hole, 1-1/2" (4 cm)




**P37486 : Keyhole Bracket**

6 mm hole, 1" (2.5 cm)




**P37487 : Pass Thru Bracket**

13 mm hole, 1" (2.5 cm)



**P37544 : Universal Keyhole Bracket**

6 positions, 7-1/2" (19 cm)



**P37545 : "L" Support Bracket**

"L" brackets to support trays

AVAILABLE

NO IMAGE  
AVAILABLE

**P37564 : Stainless Steel Riser Bracket**

Used to elevate single wide brackets from the case floor, 0.5" x 1.5"

NO IMAGE  
AVAILABLE

**P37593 : Accessory Box w/o Hardware**

Without hardware, ext. size: 1-1/2" x 3-1/2" x 1" (4 x 9 x 2.5 cm) int. size: 1.44" x 3.16" x 1" (3.5 x 8 x 2.5 cm)

NO IMAGE  
AVAILABLE

**P37872 : Lap Chole Instrument Insert**

Holds 14 instruments (fits in a 10" x 20" x 6" container), ext. size: 9.9" x 19-1/2" x 4.9"

NO IMAGE  
AVAILABLE

**P42065 : DEAN-SHALLCROSS GALL DUCT FORCEPS**

8" (20.5 cm)

NO IMAGE  
AVAILABLE

**P44924 : BAILEY FORCEPS**

Very fine angled jaws, serrated, 90 deg., 5-1/2" (14 cm)