

## Laparoscopy

NO IMAGE  
AVAILABLE

**728103 : BIPOLAR GRASPING FORCEPS**

Fenestrated, fenestrated, 5 mm, 45 cm

NO IMAGE  
AVAILABLE

**728114 : IRRIGATING/ASPIRATING NEEDLE**

15 gauge tip, Working length 13" (33cm)

NO IMAGE  
AVAILABLE

**728115 : IRRIGATING/ASPIRATING NEEDLE**

15 gauge tip, Working length 17 3/4" (45cm)

NO IMAGE  
AVAILABLE

**728116 : IRRIGATING/ASPIRATING NEEDLE**

19 gauge tip, Working length 13" (33cm)

NO IMAGE  
AVAILABLE

**728117 : IRRIGATING/ASPIRATING NEEDLE**

19 gauge tip, Working length 17 3/4" (45cm)

NO IMAGE  
AVAILABLE

**728118 : IRRIGATING/ASPIRATING NEEDLE**

22 gauge tip, Working length 13" (33cm)

NO IMAGE  
AVAILABLE

**728119 : IRRIGATING/ASPIRATING NEEDLE**

22 gauge tip, Working length 17 3/4" (45cm)

NO IMAGE  
AVAILABLE

**728120 : LAPAROSCOPIC TELESCOPE**

5mm, 30cm, 0

NO IMAGE  
AVAILABLE

**728121 : LAPAROSCOPIC TELESCOPE**

5mm, 30cm, 30

NO IMAGE  
AVAILABLE

**728122 : LAPAROSCOPIC TELESCOPE**

5mm, 30cm, 45

NO IMAGE  
AVAILABLE

**728125 : LAPAROSCOPIC TELESCOPE**

10mm, 33cm, 0

NO IMAGE  
AVAILABLE

**728126 : LAPAROSCOPIC TELESCOPE**

10mm, 33cm, 30

NO IMAGE  
AVAILABLE

**728130 : LAPAROSCOPIC TELESCOPE**

10mm, 33cm, 45

NO IMAGE  
AVAILABLE

**728131 : LAPAROSCOPIC TELESCOPE**

10mm, 45cm, 0

NO IMAGE  
AVAILABLE

**728132 : LAPAROSCOPIC TELESCOPE**

10mm, 45cm, 30

NO IMAGE  
AVAILABLE

**728133 : LAPAROSCOPIC TELESCOPE**

10mm, 45cm, 45

NO IMAGE  
AVAILABLE

**728134 : LAPAROSCOPIC TELESCOPE**

Pediatric, 2.7mm, 18cm, 0

NO IMAGE  
AVAILABLE

**728135 : LAPAROSCOPIC TELESCOPE**

Pediatric, 2.7mm, 18cm, 30

NO IMAGE  
AVAILABLE

**728140 : PALPATION PROBE**

32cm

NO IMAGE  
AVAILABLE

**728141 : PALPATION PROBE**

45cm

AVAILABLE

NO IMAGE  
AVAILABLE

**728142 : SUCTION / IRRIGATION DUAL VALVE**

5mm tip, Working length 11" (28cm)

NO IMAGE  
AVAILABLE

**728144 : SUCTION / IRRIGATION DUAL VALVE**

5mm tip, Working length 13 3/8" (34cm)

NO IMAGE  
AVAILABLE

**728145 : SUCTION / IRRIGATION DUAL VALVE**

10mm tip, Working length 13 3/4" (35cm)

NO IMAGE  
AVAILABLE

**728146 : Lazy "S" Retractor**

10 mm wide, rigid, Set of two

NO IMAGE  
AVAILABLE

**728148 : Five Finger 10mm Fan Retractor**

95 mm width blade