

**Laparoscopy**




**728149 : Five Finger 10mm Fan Articulating Retractor**




**728152 : HASSON CANNULA, 10/11MM**




**728153 : CANNULA CAP WHITE 11MM 15MM**



**728154 : CANNULA CAP BLUE 8MM 18MM**



**728155 : HASSON CANNULA, AUTO VALVE, 10/11MM**



**728157 : REDUCER SLEEVE 10/11MM TO 5MM**

NO IMAGE  
AVAILABLE

**728158 : REDUCER SLEEVE 12MM TO 5MM**

NO IMAGE  
AVAILABLE

**728159 : REDUCER SLEEVE 12MM TO 10MM**

NO IMAGE  
AVAILABLE

**728160 : REDUCER SLEEVE 5MM TO 3MM**

NO IMAGE  
AVAILABLE

**728161 : BLUNT NEEDLE ELECTRODE 33CM**

NO IMAGE  
AVAILABLE

**728162 : BLUNT NEEDLE ELECTRODE 45CM**

NO IMAGE  
AVAILABLE

**728168 : SPATULA ELECTRODE, 45CM**

NO IMAGE  
AVAILABLE

**728169 : BLUNT NOSE ELECTRODE, 33CM**



**728170 : BLUNT NOSE ELECTRODE, 45CM**



**728171 : BLUNT NEEDLE SUCTION COAGULATOR**

5mm tip, Working length 13" (33cm)




**728172 : BLUNT NEEDLE SUCTION COAGULATOR**

5mm tip, Working length 17 3/4" (45cm)



**728174 : RIGHT ANGLE (L-HOOK) SUCTION COAGULATOR**

45cm




**728176 : CURVED J-HOOK SUCTION COAGULATOR ELECTRODE**

1-3/4" (4.5 cm)



**728177 : SPATULA SUCTION COAGULATOR**

33cm



**728178 : SPATULA SUCTION COAGULATOR**

45cm

AVAILABLE

NO IMAGE  
AVAILABLE

**728179 : BLUNT NOSE SUCTION COAGULATOR**

33cm

NO IMAGE  
AVAILABLE

**728180 : BLUNT NOSE SUCTION COAGULATOR**

Working length 17 3/4" (45cm)

NO IMAGE  
AVAILABLE

**728181 : BLUNT NEEDLE IRRIGATION/SUCTION COAGULATOR**

Working length 13" (33cm)

NO IMAGE  
AVAILABLE

**728182 : BLUNT NEEDLE IRRIGATION/SUCTION COAGULATOR**

45cm

NO IMAGE  
AVAILABLE

**728184 : RIGHT ANGLE (L-HOOK) IRRIGATION/SUCTION  
COAGULATOR**

45cm