

Neurosurgery



NO IMAGE
AVAILABLE

P21504 : Yasargil Micro Scissors

Bayonet, straight, beveled micro blades, 9" (23cm)



NO IMAGE
AVAILABLE

P21512 : Yasargil Micro Scissors

Bayonet, curved blades, 9" (23cm)



NO IMAGE
AVAILABLE

P21858 : Yasargil Micro Scissors

Curved, 6 1/4" (16cm)



NO IMAGE
AVAILABLE

P22591 : Taylor Dura Scissors

5 3/4" (14.5cm)



NO IMAGE
AVAILABLE

P22593 : Dandy Trigeminal Scissors


6 3/4" (17cm)



NO IMAGE
AVAILABLE


P21499 : Yasargil Micro Scissors

Bayonet, straight blades, 8" (20.5cm)




524010 : PILLING TELESCOPE WITH FIBEROPTIC ILLUMINATION

5.5 mm, 90 lateral, 54.5 cm overall length




P14486 : SYNOVECTOMY RONGEUR

Strong curve, 5" (12.5cm)




P14619 : LISTON-MINI BONE CUTTING FORCEPS

Double-action, angled up, 6" (15 cm)




P18858 : Dandy Nerve Hook, Right

9" (23cm)




P18952 : SELLA PUNCH

90 deg. up, 1mm bite, 7" (18cm)



P18964 : SELLA PUNCH

40 deg. up, 1mm bite, 8" (20.5cm)



P19006 : KERRISON MICRO RONGEUR

90 deg. up, 9mm opening, 2mm, 7" (18cm)

NO IMAGE
AVAILABLE

P19156 : KERRISON MICRO RONGEUR

40 deg. up, 3mm, 11" (28cm)

NO IMAGE
AVAILABLE

P19158 : KERRISON MICRO RONGEUR

40 deg. up, 4mm, 11" (28cm)

NO IMAGE
AVAILABLE

P19160 : KERRISON MICRO RONGEUR

40 deg. up, 5mm, 11" (28cm)

NO IMAGE
AVAILABLE

P19162 : KERRISON MICRO RONGEUR

40 deg. up, 6mm, 11" (28cm)

NO IMAGE
AVAILABLE

P19181 : SPURLING KERRISON PUNCH

90 deg. up, 15mm opening, 3mm, 6" (15cm)

NO IMAGE
AVAILABLE

P19193 : SPURLING KERRISON PUNCH

90 deg. down, 15mm opening, 3mm, 6" (15cm)

NO IMAGE
AVAILABLE

P19279 : SPURLING KERRISON PUNCH

90 deg. down, 15mm opening, 3mm, 8" (20.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

P19458 : FERRIS-SMITH-KERRISON RONGEUR

40 deg. up, 3mm, 7" (18cm)

NO IMAGE
AVAILABLE

P19461 : FERRIS-SMITH-KERRISON RONGEUR

40 deg. up, 5mm, 7" (18cm)

NO IMAGE
AVAILABLE

P19519 : FERRIS-SMITH-KERRISON RONGEUR

40 deg. up, 2mm, 11" (28cm)

NO IMAGE
AVAILABLE

P19520 : FERRIS-SMITH-KERRISON RONGEUR

40 deg. up, 3mm, 11" (28cm)

NO IMAGE
AVAILABLE

P19576 : CUSHING IVD RONGEUR

Straight, 2 x 10mm, 9" (23cm)