

Ophthalmic




P24131 : FRANCIS FOREIGN BODY SPUD




P24171 : KELLY DESCOMET MEMBRANE PUNCH




P24180 : MEYERHOEFER CHALAZION CURETTE, #0
1.5mm



P24183 : MEYERHOEFER CHALAZION CURETTE, #3
2.5mm



P24184 : MEYERHOEFER CHALAZION CURETTE, #4
3.5mm



P24208 : Jaffe Lens Spatula
0.5mm blade

NO IMAGE
AVAILABLE

P24209 : Elsching Cyclo Spatula

1mm blade

NO IMAGE
AVAILABLE

P24211 : Castro Cyclo Spatula

Double-ended, 0.5mm wide

NO IMAGE
AVAILABLE

P24220 : Wheeler Spatula

Double-ended

NO IMAGE
AVAILABLE

P24227 : Knolle Lens Nucleus Spatula and Gauge

NO IMAGE
AVAILABLE

P24232 : BECHERT NUCLEUS ROTATOR

Angled

NO IMAGE
AVAILABLE

P24238 : DRYSDALE NUCLEUS MANIPULATOR

Paddle-shaped

NO IMAGE
AVAILABLE

P24239 : DRYSDALE NUCLEUS SPATULA

Paddle-shaped



P24254 : LESTER LENS PUSHER

Hourglass tip, angled



P24255 : LESTER IOL MANIPULATOR




P24263 : WILDER LOOP LENS

Serrated



P24270 : KIRBY EXPRESSOR HOOK AND LOOP




P24276 : BUNGE EVISCERATION SPOON

Small



P24277 : BUNGE EVISCERATION SPOON

Large



P24278 : WELLS ENUNCLEATION SPOON

Large

AVAILABLE

NO IMAGE
AVAILABLE

P24279 : FOSTER ENUCLEATION SNARE

w/ pkg. of 12 snares

NO IMAGE
AVAILABLE

P24280 : TWISTED SNARE WIRE

pkg. of 12

NO IMAGE
AVAILABLE

P24283 : Graether Micro Iris Retractor Collar Button

Angled

NO IMAGE
AVAILABLE

P24292 : Shepard Iris Hook, Small

Blunt

NO IMAGE
AVAILABLE

P24303 : Sinskey Modified Hook