

Ophthalmic



NO IMAGE
AVAILABLE

P24312 : Sinskey II Lens Manipulating Hook

Round handle, angled



NO IMAGE
AVAILABLE

P24313 : Sinskey Lens Manipulating Hook

Round handle, angled, 0.25mm dia.



NO IMAGE
AVAILABLE

P24314 : Osher Y-Hook

Round handle, straight



NO IMAGE
AVAILABLE

P24321 : Kuglen Iris Hook Lens Manipulator

Push-pull angled



NO IMAGE
AVAILABLE

P24323 : Kuglen Iris Hook Lens Manipulator

Clover leaf model



NO IMAGE
AVAILABLE

P24349 : NAGAHARA NUCLEUS CHOPPER, RIGHT HAND

For use in the right hand

NO IMAGE
AVAILABLE

P24368 : Helveston Muscle Hook

Extra-delicate, 10mm blunt point

NO IMAGE
AVAILABLE

P24370 : Vongraefe Strabismus Hook, Size #1

7mm

NO IMAGE
AVAILABLE

P24371 : Vongraefe Strabismus Hook, Size #2

9mm

NO IMAGE
AVAILABLE

P24373 : Jameson Muscle Hook, Small

9mm curve to tip, 5" (12.5cm)

NO IMAGE
AVAILABLE

P24374 : Jameson Muscle Hook, Large

NO IMAGE
AVAILABLE

P24378 : Green Strabismus Hook

NO IMAGE
AVAILABLE

P24382 : CASTRO CALIPER, 0-20MM SCALE

Straight



P24391 : 6" RULER, GRADUATED IN MM AND INCHES

Stainless steel



P24510 : Nichamin Nucleus Chopper-Spatula

Double-ended



P24694 : COLIBRI FORCEPS

Standard, 0.3 mm 1 x 2 teeth, w/ platform



P24705 : PIERSE-COLIBRI CORNEAL FORCEPS

Very delicate, curved, 0.3 mm tips, 1 x 1 teeth




P24711 : BONN SUTURING FORCEPS

0.12 mm, 1 x 2 teeth



P24721 : KIRBY CORNEAL SCLERAL FORCEPS



P24738 : CASTROVIEJO SUTURE FORCEPS

Standard handle, 0.9 mm 1 x 2 teeth

AVAILABLE

NO IMAGE
AVAILABLE

P24743 : CASTROVIEJO SUTURE FORCEPS

Wide handle, 0.9 mm 1 x 2 teeth, w/ platform

NO IMAGE
AVAILABLE

P24746 : STERN CASTROVIEJO FORCEPS

W/ thumb catch lock, 0.5 mm teeth, w/ platform

NO IMAGE
AVAILABLE

P24747 : STERN CASTROVIEJO FORCEPS

Curved, 0.5 mm teeth, w/ catch lock, w/o platform

NO IMAGE
AVAILABLE

P24757 : PAUFIQUE SUTURE FORCEPS

1 x 2 teeth, w/ platform

NO IMAGE
AVAILABLE

P24770 : CASTROVIEJO FIX FORCEPS

Wide handle, 0.5 mm teeth