

Ophthalmic




P24781 : GREEN FIX FORCEPS

5 mm wide




P24786 : BISHOP-HARMON TISSUE FORCEPS

1 x 2, 0.3 mm teeth




P24797 : DRESSING FORCEPS

Delicate, straight, non-magnetic




P24804 : TROUTMAN SUPERIOR RECTUS FORCEPS

Angled, 1 x 2 teeth



P24812 : KELMAN-MCPHER TYING FORCEPS

Standard handle, angled, smooth, 7.5 mm



P24813 : KELMAN-MCPHER TYING FORCEPS

Wide handle, angled, smooth, 7.5 mm

NO IMAGE
AVAILABLE

P24815 : SHEETS-MCPHER TYING FORCEPS

Standard handle

NO IMAGE
AVAILABLE

P24816 : BECHERT-MCPHER TYING FORCEPS

Standard handle

NO IMAGE
AVAILABLE

P24821 : UTRATA FORCEPS

Short handle

NO IMAGE
AVAILABLE

P24822 : UTRATA FORCEPS

Flat, long, wide handle

NO IMAGE
AVAILABLE

P24823 : UTRATA CAPSULORHEXIS FORCEPS

Three-hole handle

NO IMAGE
AVAILABLE

P24826 : MASKET CAPSULORHEXIS FORCEPS

11 mm shanks, 4" (10 cm)

NO IMAGE
AVAILABLE

P24836 : TENNANT TYING FORCEPS

Extra-delicate smooth jaws, straight, w/ platform



P24837 : TENNANT TYING FORCEPS

Extra-delicate smooth jaws, curved, w/ platform



P24838 : GIRARD FINE TYING FORCEPS

Straight, 4" (10 cm)




P24839 : TENNANT TYING FORCEPS

Extra-delicate smooth jaws, straight, w/ platform, knurled handle




P24852 : TYING FORCEPS

Straight, heavy shafts



P24855 : SHEPARD CURVED TYING FORCEPS

Wide handle, 4-1/4" (11 cm)



P24862 : KIRBY CAPSULE FORCEPS



P24901 : STEINERT PICKUP AND HOLDING FORCEPS

AVAILABLE

NO IMAGE
AVAILABLE

P24906 : ERNEST NUCLEUS CRACKER FORCEPS

NO IMAGE
AVAILABLE

P24910 : JAMESON MUSCLE FORCEPS, CHILD, RIGHT

NO IMAGE
AVAILABLE

P24917 : BERKE PTOSIS FORCEPS

20mm jaw, 4" (10cm)

NO IMAGE
AVAILABLE

P24918 : BERKE PTOSIS FORCEPS

27mm jaw

NO IMAGE
AVAILABLE

P24943 : SCHAAF FOREIGN BODY FORCEPS