

Ophthalmic




P25213 : O'CONNER SCLERAL MARKER




P25225 : WILDER LACRIMAL DILATOR, SIZE #2

Medium taper




P25229 : PIGTAIL PROBE

6 1/8" (15.5cm)




P25244 : WILLIAMS LACRIMAL PROBE, SIZE 00/0



P25283 : RANDOLPH CYCLO CANNULA

19g, front



P25334 : VISCOELASTIC ASPIRATING CANNULA

19 Ga.

NO IMAGE
AVAILABLE

P25335 : IRRIGATION CANNULA

27 Ga., angled tip

NO IMAGE
AVAILABLE

P25336 : AIR INJECTION CANNULA

30 Ga.

NO IMAGE
AVAILABLE

P25347 : AIR INJECTION CANNULA

27 Ga., angled 7mm from tip

NO IMAGE
AVAILABLE

P25349 : TROUTMAN ALPHA-CHYMO CANNULA

25g

NO IMAGE
AVAILABLE

P25360 : BISHOP HARMON IRRIGATION CANNULA

19 Ga.

NO IMAGE
AVAILABLE

P25377 : LACRIMAL CANNULA

Straight, 23 Ga., reinforced, 2" (5cm)

NO IMAGE
AVAILABLE

P25387 : KNOLLE-PEARCE LENS LOOP CANNULA

23 Ga.



P25394 : OSHER "J" CANNULA

27 Ga.



P25401 : WELSH CORTEX EXTRACTOR

Flat, olive tip, 25 ga., 25mm




P25478 : INFUSION HANDLE

w/ finger cut-off hole, stainless steel, 5" (12.5cm)



P25495 : CARTER SPHERE INTRODUCER




P25565 : OPEN TREPHINE HANDLE

Fits 8mm trephines



P25581 : Kuglen Iris Hook and Lens Manipulator

Angled



P25594 : NORDAN TITANIUM TYING FORCEPS

Round handle, straight, 4-3/8" (11 cm)

AVAILABLE

NO IMAGE
AVAILABLE

P25595 : NORDAN TITANIUM TYING FORCEPS

Round handle, curved, 4-1/4" (11 cm)

NO IMAGE
AVAILABLE

P25601 : TYING FORCEPS, STRAIGHT

Flat serrated handle, 7 mm platform, 4-1/2" (11.5 cm)

NO IMAGE
AVAILABLE

P37415 : CHANG PHACO CHOPPER, RIGHT HAND

For use in the right hand

NO IMAGE
AVAILABLE

P41933 : ERNEST MICRO NUCLEUS CRACKER

Extra-delicate, serrated paddle tips