

Ophthalmic



NO IMAGE
AVAILABLE

427384 : CASTROVIEJO SUTURE FORCEPS

Straight, Tying Platform, 1 x 2 Teeth, 0.12mm tip, 3 7/8" (10 cm)



NO IMAGE
AVAILABLE

423490 : Westcott Tenotomy Scissors

Flat Handle, Curved left, 4 3/4" (12.0cm)



NO IMAGE
AVAILABLE

P24355 : Fixation Hook, Small

Double sharp, 120mm



NO IMAGE
AVAILABLE

P25166 : Jaegar Lid Plate Retractor

Stainless steel



NO IMAGE
AVAILABLE

P24380 : Gass Retinal Detachment Hook

5 5/8" (14.5cm)



NO IMAGE
AVAILABLE

P25178 : Schepens Orbital Retractor

5 7/8" (15cm)

NO IMAGE
AVAILABLE

P24199 : Culler Iris Spatula

2 mm wide

NO IMAGE
AVAILABLE

425590 : COLIBRI CORNEAL TISSUE FORCEPS

1 x 2 Teeth, 0.12mm tip, 3" (8 cm)

NO IMAGE
AVAILABLE

P24300 : Sinskey Lens Manipulating Hook

Straight, blunt, 0.25mm dia.

NO IMAGE
AVAILABLE

421640 : Desmarre Lid Retractor, Size 4

17mm width blade, 5 3/8" (13.5cm)

NO IMAGE
AVAILABLE

P25172 : Jaffe Lid Retractor, Large

Pair

NO IMAGE
AVAILABLE

P24841 : TENNANT TYING FORCEPS

Extra-delicate smooth jaws, curved, w/ platform, knurled handle

NO IMAGE
AVAILABLE

144340 : Knapp Iris Scissors

Curved, Sharp/Blunt blades, 4 1/2" (11.5cm)



P25013 : CASTROVIEJO NEEDLE HOLDER

Delicate, straight, 10mm jaws



P25015 : CASTROVIEJO NEEDLE HOLDER

Delicate, curved, 10mm jaws, w/lock



P25017 : CASTROVIEJO NEEDLE HOLDER

Straight, standard jaws, w/lock



P25018 : CASTROVIEJO NEEDLE HOLDER

Curved, standard jaws, w/lock




P24534 : Troutman-Castroviejo Corneal Scissors, Mini, Left



P25023 : CASTROVIEJO NEEDLE HOLDER

Delicate, curved, 9mm jaws, w/lock



P25034 : COHAN MICRO NEEDLE HOLDER

Curved, w/ lock

AVAILABLE

NO IMAGE
AVAILABLE

P25035 : COHAN MICRO NEEDLE HOLDER

Curved, w/o lock

NO IMAGE
AVAILABLE

P25052 : BARRAQUER NEEDLE HOLDER

Delicate, curved, w/o lock, 5" (12.5cm)

NO IMAGE
AVAILABLE

P25058 : BARRAQUER NEEDLE HOLDER

Extra delicate, straight, w/o lock, 5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

425200 : MCPHERSON TYING FORCEPS

Straight, Smooth, 3 3/8" (9 cm)

NO IMAGE
AVAILABLE

P25060 : BARRAQUER NEEDLE HOLDER

Extra delicate, curved, w/o lock, 5 1/4" (13.5cm)