


Orthopedic

**KM47904 : SMALL QUICK COUPLING HANDLE**

for small and large screw sets, 4 1/2" (11.5cm)

**KM48079 : DRILL BIT, JACOBS CHUCK**

3.2mm dia., 70mm , 7 1/8" (18cm)

**KM35292 : Super-Cut Mayo Scissors**

Straight, 6 3/4" (17.0cm)

**KM48054 : DRILL BIT QUICK COUPLING**

2.0mm dia., 50mm, 3 7/8" (10cm), 75mm working length

**594402 : HOEN HOOK**

Long prong, 9-1/4" (23.5 cm)

**KM46906 : Fukuda Style Shoulder Retractor, Large**

38mm width x 81mm length blade, 7" (18.0cm)



KM47880 : SHARP HOOK

6" (15.2 cm)



KM48240 : MICHELE TREPHINE

7 1/2" (19cm), 4mm ID / 6mm OD



KM47936 : SMALL SCREW DEPTH GAUGE

6 1/4" (16cm), up to 50mm



KM51162 : Meyering Laminectomy Retractor, Medium

sharp teeth, 25mm width blade, 9 1/2" (24.0cm)



P30366 : CORKSCREW T-HANDLE

7 1/2" (19cm)



KM35370 : Super-Cut Metzenbaum Scissors

Curved, 4 1/2" (11.5cm)



KM47902 : Mini Quick Coupling Handle

For holding quick coupling drills, countersinks, taps and screwdrivers. For mini screw set.
4-1/8" (10.5 cm)

KM49634 : KEY ONLY FOR KM49630

3 3/8" (8.5cm)



0650110 : HUDSON BURR

4" (10.0cm), 9mm



KM47906 : BALL SPIKE PELVIC STRAIGHT

11 3/4" (30.0cm)




KM48342 : BONE HOOK, SMALL

Sharp, 6 mm x 10 mm hook, 9" (22.9 cm)



KM35490 : Super-Cut Stevens Scissors

Curved, 4 1/2" (11.5cm)




KM46602 : PERCY AMPUTATION SHIELD

12" (30.5cm)



KM46852 : Hohmann Retractor

22mm width blade, 16" (40.5cm)



NO IMAGE

AVAILABLE



KM48158 : FRAZIER HOOK

Sharp, curved, 2.5 mm diam., single prong, 5" (12.7 cm)



KM52765 : GELPI-SELETZ RETRACTOR

6 1/4" (16cm)




KM48058 : DRILL BIT QUICK COUPLING

2.7mm dia., 29mm, 3 7/8" (10cm), 75mm working length



KM48226 : LAMINA SPREADER FLAT BLADES

10 3/4" (27.4cm)



KM35394 : Super-Cut Metzenbaum Scissors

Curved, 7" (18.0cm)