

Orthopedic



028101 : LEMPERT RONGEUR

Curved, delicate, 10 x 3mm




051082 : LAMBOTTE OSTEOTOME

Straight, 2mm tip width, 5" (12.5cm)




051084 : LAMBOTTE OSTEOTOME

Straight, 4mm tip width, 5" (12.5cm)




051086 : LAMBOTTE OSTEOTOME

Straight, 6mm tip width, 5" (12.5cm)



051088 : LAMBOTTE OSTEOTOME

Straight, 8mm tip width, 5" (12.5cm)



051090 : LAMBOTTE OSTEOTOME

Straight, 10mm tip width, 5" (12.5cm)

NO IMAGE
AVAILABLE

051092 : LAMBOTTE OSTEOTOME

Straight, 12mm tip width, 5" (12.5cm)

NO IMAGE
AVAILABLE

051095 : LAMBOTTE OSTEOTOME

Straight, 15mm tip width, 5" (12.5cm)

NO IMAGE
AVAILABLE

051102 : Hex-Handle Osteotome

2mm tip width, flat, thin blade. 6" (15 cm).

NO IMAGE
AVAILABLE

051103 : Hex-Handle Osteotome

3mm tip width, flat, thin blade. 6" (15 cm)

NO IMAGE
AVAILABLE

051104 : Hex-Handle Osteotome

4mm tip width, flat, thin blade. 6" (15 cm)

NO IMAGE
AVAILABLE

051105 : Hex-Handle Osteotome

5mm tip width, flat, thin blade. 6" (15 cm)

NO IMAGE
AVAILABLE

051106 : Hex-Handle Osteotome

6mm tip width, flat, thin blade. 6" (15 cm)



051107 : Hex-Handle Osteotome

7mm tip width, flat, thin blade. 6" (15 cm)



051109 : Hex-Handle Osteotome

9mm tip width, flat, thin blade. 6" (15 cm)




051110 : Hex-Handle Osteotome

10mm tip width, flat, thin blade. 6" (15 cm)




051150 : COTTLE OSTEOTOME

Straight, 4mm tip width, 7 1/4" (18.5cm)




051152 : COTTLE OSTEOTOME

Straight, 6mm tip width, 7 1/4" (18.5cm)



051155 : COTTLE OSTEOTOME

Straight, 9mm tip width, 7 1/4" (18.5cm)



051158 : COTTLE OSTEOTOME

Straight, 12mm tip width, 7 1/4" (18.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

051190 : ANDERSON-NIEVERT OSTEOTOME

Curved right, 8" (20.5cm)

NO IMAGE
AVAILABLE

051191 : ANDERSON-NIEVERT OSTEOTOME

Curved left, 8" (20.5cm)

NO IMAGE
AVAILABLE

051195 : ANDERSON-NIEVERT OSTEOTOME

Straight, 8" (20.5cm)

NO IMAGE
AVAILABLE

051500 : Mallet LEAD FILLED HEAD

7"

NO IMAGE
AVAILABLE

051510 : SOLID COPPER HEAD Mallet

Solid, Head diameter 3/4" (20mm), 6 1/2" (16.5cm)